
ULINE H-5101

SHOWER EYEWASH STATION

1-800-295-5510
uline.com

INSTALLATION

- Assemble unit in accordance with the manufacturer's drawing.
- Install in close proximity to hazard and in a location clearly identified and without obstructions.
- Place shower head 82-96" from the surface at which the user stands.
- Position center of spray pattern at least 16" from any obstruction.
- Arrange spray heads 33-45" from the surface at which the user stands and a minimum of 6" from any obstruction.
- Connect unit with a recommended minimum 1¼" NPT line to a tepid, potable water source capable of supplying 30-80 psi flowing pressure for a minimum of 15 minutes.
- If shut valves are installed, inline provisions should be made to prevent unauthorized shutoff.
- Protect unit from freezing.
- Check for any leaking while all valves are in the closed position.
- Keep valves in the open position upon activation.
- With the shower valve in the fully opened position, spray pattern should be a minimum of 20" in diameter at 60" above the surface at which the user stands. Determine flow rate is at least 20 GPM and the water is tepid.
- With eyewash valve in fully opened position and using a test gauge, the spray pattern should meet required pattern. Determine flow rate is at least 0.4 GPM and the water is tepid.

ASSEMBLY

1. Assemble the base plate (A), 1¼" IPS x 6¼" nipple (B) and plugged outlet tee (C).

IMPORTANT! Plugged outlet tee (C) must be oriented with plugged end down to ensure proper drainage.

2. Assemble 1¼" IPS x 16½" nipple (D) and frame assembly (E). Position frame assembly for proper access by users.
3. Assemble 1¼" IPS x 32" nipple (F) and inlet tee (G). Position tee for connection to supply. If the top vertical inlet will be employed, plug tee (G) with supplied 1¼" IPS square head plug.
4. Assemble 1¼" IPS x 22" nipple (H) and 1¼ x 1" IPS reducing inlet tee (I). Position tee forward for proper location of shower. If the back horizontal inlet will be employed, plug tee (I) with supplied 1¼" IPS square head plug.
5. Assemble 1¼" IPS x 8⅞" nipple (J), emergency shower valve (K), 1" IPS x 10½" nipple (L), elbow (M) and shower head (N).
6. Place cupped washer (O), rubber gasket (P) and bowl (Q) onto waste receptor. Center components to ensure a watertight seal.
7. Thread drain plate (R) down onto mounting shank as far as possible. If necessary, rotate bowl clockwise to draw drain plate down onto mounting shank.
8. Thread spray head assembly (S) onto mounting shank. Position outlet heads for proper access by user.
9. Secure unit to the floor using the four mounting slots in base plate (A) and floor mounting hardware (not included).

NOTE: Plumbing is required for the water supply.

10. Attach orange plastic paper wrap by prying slit open and wrapping it on the pipe.

NOTE: All threads are tapered pipe and require pipe joint compound or thread sealing tape to ensure watertight connection.

NOTE: Test the entire unit for proper operation.

NOTE: The use of dielectric unions is required when connecting to copper supply to avoid premature corrosion.

NOTE: Pull rod to start shower. Push pedal to start eyewash.

Figure 1




Figure 2

OPERATION

SHOWER

1. To activate, pull down on the triangular pull rod in one swift motion. The pull rod will open the 1" ball valve and allow water to flow through the shower head. (See Figure 3)


NOTE: The valve will stay open and allow the shower to operate "hands-free" and continuously until the valve is manually closed.

2. To deactivate, push up on the pull rod to return it to the stored position. This will close the 1" ball valve and turn off the shower.

EYEWASH

1. To activate, push the pedal back in one swift motion. The pedal will open the 1/2" ball valve and allow water to flow through the eyewash. (See Figure 4)

NOTE: The valve will stay open and allow the eyewash to operate "hands-free" and continuously until the valve is manually closed.

2. To deactivate the eyewash, pull the pedal back to the original position. This will close the 1/2" ball valve and turn off the eyewash.

MAINTENANCE

WEEKLY ACTIVATION

IMPORTANT! ANSI Z358.1-2014 requires that emergency equipment be activated at least weekly to ensure proper operation. Weekly inspections are required to verify performance standards listed below:

- Ensure that there is sufficient water delivered to the emergency unit for effective operation.
- Flush out any buildup of sediment or other impurities in the water lines that could prevent the unit from operating effectively.


- Minimize microbial contamination due to stagnant water standing in the piping and the unit.
- Verify that the equipment is free of obstructions.
- Ensure that the unit delivers adequate flow.

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