

**ULINE** H-5276, H-5278  
 H-5279, H-6097  
**OHAUS ADVENTURER®**  
**BALANCE SCALE**

1-800-295-5510  
 uline.com

**OVERVIEW OF CONTROLS**

**CONTROL PANEL**

**CONTROL PANEL PARTS**

#	DESCRIPTION
1	ON/OFF Button
2	Print Button
3	Zero Button
4	Calibration Button
5	Tare Function Buttons

**CONTROL FUNCTIONS**

BUTTON					
<b>PRIMARY FUNCTION (Short Press)</b>	<b>ON/OFF</b> Turns scale on if scale is off.	<b>PRINT</b> Prints the present data to a computer printer.	<b>ZERO</b> Performs zero operation.	<b>CALIBRATION</b> Performs calibration operation.	<b>TARE</b> Performs tare operation.
<b>SECONDARY FUNCTION (Long Press)</b>	If scale is on, turns scale off.  <b>NOTE: The balance will automatically power on when power is connected.</b>				

**MAIN APPLICATION SCREEN**

Application — Dynamic Weighing

Instructional Messages — Place sample on the pan. Press Start to begin averaging.

Stability (\*), Net (NET), Gross (G) and/or center of zero (>0<) indicators — \* >0<

Reference Fields — Gross: 0.00 g Net: 0.00 g Averaging Time: 5 s Automatic Mode: Off

Capacity and Readability — Max 4201 g d=0.01 g e=0.01 g

Result Field: Information varies by application. Touch g to change unit. — 0.00 g

Application Buttons: Functions vary by application. — Item Settings Start Save to USB Menu

# OPERATION

## TURNING BALANCE ON/OFF


- To turn the balance ON, press the  button.
- To turn the balance OFF, press and hold the  button.

## MENU NAVIGATION

### TO ENTER/EXIT THE MENU


All menu navigation is performed by touching the display. To enter Menu, touch **MENU** from any application home screen. The main menu appears, with buttons for **BACK** and **EXIT**. (See Figure 1)

Figure 1

Continue touching the appropriate list item to navigate to the menu items. (See Figure 2)

Figure 2


## CHANGING MENU SETTINGS

To change a menu setting, navigate to that setting using the following steps:

- Enter the Menu:** From any application screen, touch **MENU**. The main menu list appears on the display. (See Figure 3)
- Select the Sub-Menu:** Find the item of the main menu list and touch it. The sub-menu appears.

Figure 3

- Select the Menu Item:** Continue until the desired setting is chosen in the menu list. Touch the setting to change it. The changed setting will be displayed as highlighted yellow for about one second to confirm the changed value.
- Exit the Menu and Return to Current Application:** After the setting is confirmed, touch **EXIT** to return to the Application.

**NOTE:** At any time the **BACK** and **EXIT** buttons can be touched to navigate to the desired area of the menu or return to the current Application. Continue until the desired setting is selected in the menu list.

## MENU LIST

Touch **MENU** to open the menu list.

**Calibration:**  
Touch to view calibration options.

**Balance Setup:**  
Touch to view and change balance settings.

**Weighing Units:**  
Touch to view and change weighing units.

**Data Maintenance:**  
Touch to view data maintenance settings.

**Communication:**  
Touch to view COM device settings and print settings.

**GLP and GMP Data:**  
Insert user data for traceability.

**Factory Reset:**  
Touch to do a factory reset of menu settings.

**Lockout:**  
Touch to view lockout options.

# OPERATION CONTINUED

## ENABLING UNITS OF MEASURE

Enter the weighing units sub-menu in the main menu to activate the desired units of measure. This menu can also be accessed by pressing the unit symbol in an application home screen.

**NOTE:** Due to national laws, the balance may not include some of the units of measure listed.

### UNITS SUB-MENU

<b>mg</b> Milligram	<b>g</b> Gram	<b>kg</b> Kilogram	<b>ct</b> Carat	<b>OZ</b> Ounce	<b>ozt</b> Troy Ounce
<b>lb</b> Pound	<b>dwt</b> Pennyweight	<b>Grain</b> Grain	<b>N</b> Newton	<b>mom</b> Momme	<b>msg</b> Mesghal
<b>Hkt</b> Tael (HK)	<b>SGt</b> Tael (SG)	<b>TWt</b> Tael (TW)	<b>tical</b> Tical	<b>tola</b> Tola	<b>bht</b> Baht
<b>Cl</b> Custom Unit 1					

**NOTE:** Only H-5276 converts to mg.

**NOTE:** If Approved mode is on, some units will not be displayed.

Use the Custom Unit to display weight in an alternative unit of measure.

The custom unit is defined using a conversion factor, where the conversion factor is the number of custom units per gram expressed in scientific notation (Factor x 10<sup>Exponent</sup>). The custom unit's name can be up to three characters.

**For example:** To display weight in troy ounces (0.03215075 troy ounces per gram), enter a factor of 0.3215075 and an exponent of -1.

## CALIBRATION

**NOTE:** The calibration unit is always in grams.

**CAUTION!** Do not disturb the balance during any calibration.

### CALIBRATION SUB-MENU

 Span Calibration	 Linearity Calibration	 Cal Test
--	---	--

**NOTE:** H-5276, H-5278, H-5279 and H-6097 have Span Calibration, Linearity Calibration and Cal Test.

## INITIAL CALIBRATION

When the balance is first installed, and when it is moved to another location, it must be calibrated to ensure accurate weighing results. The balance can be manually calibrated with external masses. Have the appropriate calibration masses available before beginning calibration.

CALIBRATION WEIGHTS			
ULINE MODEL NO.	OHAUS MODEL NO.	SPAN CAL POINTS (G)	LINEARITY CAL POINTS (G)
H-5276	AX423/E	0, 100, 200, 300, 400	0, 200, 400
H-5278	AX1502/E	0, 500, 1,000, 1,500	0, 1,000, 1,500
H-5279	AX8201/E	0, 2,000, 4,000, 6,000, 8,000	0, 4,000, 8,000
H-6079	AX822/E	0, 200, 400, 800	0, 400, 800

## SPAN CALIBRATION (SINGLE POINT)

Span calibration uses two calibration points, one at zero load and the other at specified full load (span).

 **NOTE:** Factory default settings are shown in bold.

With the balance turned on and no load on the pan, touch **SPAN CALIBRATION** to initiate the procedure. Additional calibration values to be used are shown on the display. The best accuracy is achieved using the mass closest to the full span value.

 **NOTE:** To change the span calibration point, touch the alternate weight shown on the display. Follow the screen instructions and place the specified calibration weight on the scale when prompted to do so. When complete, the display shows the span calibration status and returns to the current application.

## LINEARITY CALIBRATION (MULTI-POINT)

Linearity calibration uses three calibration points, one at zero load and the others at specified loads. With no load on the scale, press **LINEARITY CALIBRATION** to begin the process. The balance captures the zero point, and then prompts for the next weight. Continue to follow the instructions on the display until the calibration is completed.

 **NOTE:** To cancel at any time, press **CANCEL**.

## CALIBRATION TEST

Use Calibration Test to compare a known calibration weight against the stored span calibration data.

With no load on the scale, press **CAL TEST** to begin the process. The balance captures the zero point, then prompts for the span value.

The display shows status, followed by the difference between the current calibration weight and the stored calibration data.

## ZERO

Press **ZERO** to zero the balance.

## TARING

With no load on the pan, press **ZERO** to set the display to zero. Place an empty container on the pan and press **TARE**. Add material to the container and its net weight is displayed. Remove container and container's weight appears as a negative number. Press **TARE** to clear.

## APPLICATION MODES

Touch the top left application field (weighing in the example below) to select the application mode. (See Figures 4-5)

Figure 4



Figure 5

**NOTE:** Before using any application, ensure the balance has been leveled and calibrated

## WEIGHING

Use this application to determine the weight of items in the selected unit of measure.

1. In the upper left portion of the home screen, (See Figure 6) select **WEIGHING** (this application is the default).
2. Press **TARE** or **ZERO** if necessary to begin.
3. Place objects on the pan to display the weight. When stable, the \* appears.
4. The resulting value is displayed in the main display line in the active unit of measure.



Figure 6

## ITEM SETTINGS

To view or adjust the current settings, touch the **ITEM SETTINGS** button. The Settings screen appears. (See Figure 7)

Figure 7


**Minimum Weight:** Establish a minimum weight value to be used to verify a reading. If an actual weight is below the established minimum weight value, it is flagged by a color change: **yellow**.

To adjust the minimum weight value, touch the **MINIMUM WEIGHT** button.


A numeric input window appears. (See Figure 8) Use the keys to enter the desired minimum weight, then press **SAVE**. The display reverts to the previous screen. To return to the Weighing home screen, touch **EXIT** at the bottom of that screen.

Figure 8


**Capacity Bar:** When set to "on", a capacity bar is displayed in the reference field. The capacity guide will show the current weight as a percentage of balance capacity. (See Figure 9)

Figure 9

If capacity bar is set to "off", the reference field will show minimum weight and sample name. (See Figure 10)

Figure 10

**Weighing Units:** Change the displayed unit.

**NOTE:** Touching the weighing unit from application home screen will also open the Weighing Units screen.

**Filter Level:** Change filtering level.

**GLP & GMP Data:** Change Good Laboratory Practices (GLP) and Good Manufacturing Practices (GMP) data.

**Print Settings:** Change printing settings.

**Level Assist:** Instructions on how to move the balance feet to level the balance.

## SAMPLE NAME

Press this button to add a sample name. An alphanumeric input window appears. Press to alternate between lower and upper case characters. (See Figures 11-12)

Figure 11 Lower Case:



Figure 12 Upper Case:


Key in the desired sample name and press **SAVE** to save the name and return to weighing home screen.

# OPERATION CONTINUED

## SAVE TO USB

Insert the USB flash drive into the USB slot located on the front of the balance. Next, press the **SAVE TO USB** button to save the data to the USB flash drive. Once saved, the button will momentarily change color to orange. (See Figure 13)

**NOTE:** The first time a USB flash drive is connected to the balance, there might be some delay time before the **SAVE TO USB** button works. This is due to the balance creating the necessary folders on the USB flash drive where the data will be stored.

The density determination and check weighing applications do not have a **SAVE TO USB** button.

**CAUTION!** The weighing data will be saved to USB every day. However, if different weighing modes are used, the data will be separately saved to individual files.

Depending on the USB drive used, all data might not be transferred from the balance or the display might freeze. If this happens, unplug the USB flash drive and try another USB flash drive.

Uline takes no responsibility if data on USB flash drive is erased or if the USB flash drive breaks while it is connected to the balance. To minimize the risk of problems arising, Uline suggests using a high quality USB flash drive.

## PARTS COUNTING

Use this application to count samples of uniform weight.

1. In the upper left portion of the home screen, select Parts Counting.
2. Press **TARE** or **ZERO** if necessary to begin.
3. Place objects on the pan to display the weight. When stable, the \* appears.

4. The resulting value is displayed in the main display line in pieces (PCS). (See Figure 14)



Figure 14


## ITEM SETTINGS

To view or adjust current settings, touch the **ITEM SETTINGS** button. The settings screen appears. (See Figure 15)

**Samples:** The sample size can be 1 to 10,000 pieces. The default sample size is 10. Once a sample size is changed, the balance will immediately open the recalculate APW screen, expecting to establish a new APW.

1. To adjust the sample size, touch the **SAMPLES** button. A numeric input window appears. (See Figure 16)


2. Use the keys to enter the desired sample size, then press **SAVE**.

# OPERATION CONTINUED


- The next screen appears, with the message to place the reference weight on the pan. (See Figure 17)

Figure 17


**Set a Known Average Piece Weight (APW):** To adjust the APW value directly, touch the **APW** button. A numeric input window appears. (See Figure 20)

Figure 20

- Place the reference weight on the pan, then touch **ACCEPT** to capture the value; the screen shows number of pcs. (See Figure 18)

Figure 18


Key in the desired APW Weight, then press **SAVE**.

The display returns to the application home screen with the new APW value in the reference field. (See Figure 19)

**Set a New Average Piece Weight (APW) – Derived:** To establish a new APW, place the reference weight on the pan and press **ACCEPT** to continue. (See Figure 21)

Figure 21

**Establish an Average Piece Weight (APW):** Each time a new type of part is counted, the nominal weight of one piece (APW) must be established using a small quantity of pieces. This APW is stored until replaced by another APW.

There are two methods to establish the APW value:

- The actual APW is known.
- The APW must be derived by weight. For this case the current sample size will be used. (See Figure 19)

 **NOTE:** The sample size that is displayed will be used. To use a different sample size, change that first. (See Figure 21)

The home screen shows the number of pieces at the new APW. (See Figure 22)

Figure 19



Figure 22

# OPERATION CONTINUED

**Auto Optimization:** When set to *ON*, improves counting accuracy by recalculating the piece weight automatically as parts are added. Default is *OFF*.

 **NOTE:** APW optimization occurs only when the number of pieces added to the pan is between one and three times the number already on the pan. If the APW is entered manually by the numeric keypad, APW auto optimization will not occur.

**Print Settings:** Change printing settings.

 **NOTE:** The *SAVE TO USB* button will only appear after a USB flash drive has been connected to the balance. See page 6 for more information.

## PERCENT WEIGHING

Use Percent Weighing to measure the weight of a sample displayed as a percentage of a preestablished reference weight. The default (or last) reference weight is displayed.

1. In the upper left portion of the home screen, select Percent Weighing. (See Figure 23)
2. Place an object on the pan. The difference between the sample and the reference weight is displayed as a percentage.



Figure 23

## ITEM SETTINGS

To view or adjust the current settings, touch the **ITEM SETTINGS** button. The Settings screen appears. (See Figure 24)


Figure 24

**Recalculate Reference Weight:** To establish a new reference weight value, touch the **RECALCULATE REFERENCE WEIGHT** button.

Follow the screen instructions to establish a new reference weight. (See Figure 25)

Figure 25

Alternatively, press the **REF. WEIGHT** button from the Percent Weighing Recalculate Ref. Wgt screen to establish a new reference weight manually through a numerical keypad.

**Print Settings:** Change printing settings.

 **NOTE:** The *SAVE TO USB* button will only appear after a USB flash drive has been connected to the balance. See page 6 for more information.

# OPERATION CONTINUED

## DYNAMIC WEIGHING

Use this application to weigh an unstable load, such as a moving animal. Two different start/reset modes can be selected: **Manual** (start and stop via key press) and **Automatic** (start and stop automatically).

### DYNAMIC WEIGHING – MANUAL (DEFAULT)

1. In the upper left portion of the home screen, select Dynamic Weighing. (See Figure 26)
2. Place objects on the pan and press the **START** button.


Figure 26


3. The balance begins a countdown (averaging process). During the countdown, the information line displays the time remaining. If necessary, press **STOP** to quit.
4. When the countdown ends, the result line is displayed and held. Press **CLEAR** to clear the held weight and return to the initial screen.

 **NOTE:** Clear the pan before beginning a new dynamic weighing cycle.

### ITEM SETTINGS

To view or adjust the current settings, touch the **ITEM SETTINGS** button. The Settings screen appears. (See Figure 27)


Figure 27

**Averaging Time:** Set the averaging time to a value between 1 and 99 seconds. Default is 5 seconds.

**Automatic Mode:** When set to ON, the cycle begins when an object is placed on the pan, and the held value is automatically reset when the object is removed from the pan. (See Figure 28)

Figure 28

**Sample Name:** Assign a name to the sample.

**Print Settings:** Change printing settings.

 **NOTE:** The **SAVE TO USB** button will only appear after a USB flash drive has been connected to the balance. See page 6 for more information.

## DENSITY DETERMINATION

Use this application to determine an object's density. (See Figure 29)

Four types of density determination can be made:

1. Solids more dense than water
2. Solid less dense than water
3. Liquid density
4. Porous material (impregnated with oil)



Figure 29

# OPERATION CONTINUED

## MEASURING THE DENSITY OF A SOLID USING WATER (DEFAULT)

Press the **ITEM SETTINGS** button to open the Density Determination settings screen. (See Figure 30)

Figure 30

Confirm the following Setups are selected:


- Density Type: Solid
- Auxiliary Liquid: Water
- Porous Material: Off

1. To adjust the water temperature value, touch the **WATER TEMP.** button.

**NOTE:** The balance calculates water density based on the water temperature value entered.

2. Measure the actual water temperature using a precision thermometer. A numeric input window appears. (See Figure 31)


Figure 31

3. Enter the actual water temperature, then press **SAVE**.


4. To return to the Density Determination home screen, touch **BACK**. (See Figure 32)

Figure 32


**Step 1 – Weigh the Sample in Air:** Press **START**. Follow screen instructions, then press **ACCEPT** to store the dry sample weight (in air). (See Figure 33)

Figure 33


**Step 2 – Weigh the Sample Submerged in Liquid:** Follow the screen instructions, then press **ACCEPT** to store the sample weight (submerged in liquid). (See Figure 34)

Figure 34

Once the necessary weights have been determined, the density of the sample is displayed in g/cm<sup>3</sup> (along with the weight in air, weight in water) on the application screen. (See Figure 35)

Figure 35


5. Press **START** to reset all the weight values and restart the process.

# OPERATION CONTINUED

## MEASURING THE DENSITY OF A BUOYANT MATERIAL USING WATER

Press the **ITEM SETTINGS** button to open settings screen. (See Figure 36)

Figure 36

Confirm the following **Setups** are selected:

- **Density Type: Solid**
- **Auxiliary Liquid: Water**
- **Porous Material: Off**

Press **BACK** to return to the density determination home screen.

Follow the same procedure as solid material above, except in density determination step two, push the sample down into the liquid until it is fully submerged.


## MEASURING THE DENSITY OF A SOLID USING AN AUXILIARY LIQUID

To enable this feature, enter the density determination setup menu and select the following:

- **Density Type: Solid**
- **Auxiliary Liquid Type: Other**
- **Porous Material: Off**

Confirm the default values displayed (liquid density, etc) are correct. (See Figure 37)


Figure 37

1. To adjust the liquid density value, touch the **AUXILIARY LIQUID** button and select **Other**.

2. A numeric input window appears. (See Figure 38)

Figure 38

3. Key in the density in **g/cm<sup>3</sup>**, then press **SAVE**. The display returns to the previous screen. To return to the density determination home screen, touch **BACK**.
4. Begin the density determination process as per above.

## MEASURING THE DENSITY OF A LIQUID USING A CALIBRATED SINKER (not supplied)

To enable this feature, enter the density setup menu and select the following; **Density Type: Liquid**. (See Figure 39)

**NOTE:** When the density type is set to **Liquid**, the liquid type and porous material selections are disabled.



Figure 39

Confirm that the default value displayed (sinker volume) is correct. To edit the default values, touch **SINKER VOLUME**.

1. To adjust the sinker volume value, touch the **SINKER VOLUME** button. (See Figure 40)


Figure 40

# OPERATION CONTINUED

2. A numeric input window is displayed. (See Figure 41)


Figure 41

3. Key in the desired value, then press **SAVE**. The display returns to the previous screen with the new value highlighted.


**NOTE:** To return to the density determination home screen, touch **BACK**.

4. Press **START** to start the process.

**Step 1 – Weigh the Sinkers in Air:** Follow the screen instructions, then press **ACCEPT** to store the sinker weight (in air). (See Figure 42)

Figure 42

**Step 2 – Weigh the Sinkers Submerged in the Test Liquid:** Follow the screen instructions, then press **ACCEPT** to store the sinker weight (submerged in liquid). (See Figure 43)


Figure 43

Once the necessary weights have been determined, the density of the liquid sample is displayed in  $\text{g/cm}^3$  (along with the weight in air, weight in water) on the application screen.

5. Press **START** to reset all the weight values and restart the process. (See Figure 44)

Figure 44

## MEASURING THE DENSITY OF POROUS MATERIAL USING OIL

To enable this feature, enter the density determination Setup menu (See Figures 45-46), and set the following:

- Density Type: Solid
- Auxiliary Liquid Type: Water
- Porous Material: On

Figure 45



Figure 46


The **DENSITY DETERMINATION – POROUS** home screen  
 Main Display Line  
 Reference Fields  
 Function Buttons

# OPERATION CONTINUED

Confirm the default values displayed (Water Temp.) are correct.

1. To edit the default values, touch **ITEM SETTINGS**. The settings screen appears. (See Figure 47)

Figure 47

The balance calculates water density based on the water temperature value entered (look-up table). Measure the actual water temperature using a precision thermometer.

2. To adjust the water temperature or oil density values, touch the **WATER TEMP.** or **OIL DENSITY** button. Numeric input windows appear. Key in the desired value, then press **SAVE**. (See Figures 48-49)

Figure 48



Figure 49

3. The display returns to the previous screen with the new value highlighted. To return to the density determination home screen, touch **EXIT**.
4. Press **BEGIN DENSITY CALCULATION**.

## STEP 1 – WEIGH THE UN-OILED SAMPLE IN AIR.

Follow the screen instructions then press **ACCEPT** to store the dry sample weight (in air). (See Figure 50)

Figure 50

## STEP 2 – WEIGH THE OILED SAMPLE IN AIR.

Follow the screen instructions then press **ACCEPT** to store the sample weight (oiled). (See Figure 51)


Figure 51

## STEP 3 – WEIGH THE OILED SAMPLE SUBMERGED IN LIQUID.

Follow the screen instructions then press **ACCEPT** to store the oiled sample weight (submerged in liquid). (See Figure 52)

Figure 52


# OPERATION CONTINUED

Once the necessary weights have been determined, the density of the sample is displayed in g/cm<sup>3</sup> (along with the weight in air, un-oiled and oiled, and weight in water) on the application screen.

The value stays on the display until **START** is touched. Press **START** to reset all the weight values and restart the process. (See Figure 53)

Figure 53

## CHECK WEIGHING

Check Weighing is used to compare the weight of a sample against target limits.

1. In the upper left portion of the home screen, select Check Weighing. (See Figure 54)
2. The default (or last) check weight limits are displayed.
3. Place objects on the pan.
4. The Under/Accept/Over status is shown in the progress bar area while the actual weight of the item is shown on the main display line.


Figure 54

To set the over limit value, touch the **OVER LIMIT** button. To set the under limit value, touch the **UNDER LIMIT** button. A numeric input window appears. (See Figures 55-56)

Figure 55



Figure 56


Enter the desired limit weight, then press **SAVE**. To return to the Check Weighing home screen, touch **EXIT**.

## ITEM SETTINGS

To view or adjust the current settings, touch the **ITEM SETTINGS** button. The settings screen appears. (See Figure 57)

Figure 57

**Sample Name:** Assign a name to the sample.

**Print Settings:** Change printing settings.

## DISPLAY HOLD

Two modes are available: (See Figure 58)

**Display Hold** – allows the user to capture and store a stable weight.

**Peak Hold** – allows the user to capture and store the highest stable weight.



Figure 58

## DISPLAY HOLD

1. In the upper left portion of the home screen, select Display Hold. (See Figure 59)
2. Place the sample on the pan and press **HOLD** at any time while the weight is being captured.
3. The main display line now shows the first stable weight.
4. Press **CLEAR** to remove the hold and return to Display Hold home screen.


Figure 59

## PEAK HOLD

1. In the upper left portion of the home screen, select Display Hold. (See Figure 60)
2. Choose Peak Hold mode in item settings. (See below)
3. Place sample on the pan and press **START** to begin.
4. Continue to weigh samples. The highest stable weight will be held.
5. To remove the hold and return to normal operation press **STOP**.


Figure 60

## ITEM SETTINGS

To view or adjust the current settings, touch the **ITEM SETTINGS** button. The settings screen appears. (See Figure 61)

Figure 61

**Mode:** Choose between Peak Hold and Display Hold (default).

**Sample Name:** Assign a name to the sample.

**Print Settings:** Change printing settings.

**NOTE:** The **SAVE TO USB** button will only appear after a USB flash drive has been connected to the balance.

## TOTALIZATION

Totalization measures the cumulative weight of a sequence of items. (See Figure 62) The cumulative total may exceed the capacity of the balance. The maximum number of samples (n) is 99.


Figure 62

1. In the upper left portion of the home screen, select **Totalization**.
2. Place item on the pan to begin. The sample weight is shown on the main display line.
3. Press **ACCUMULATE** to add the weight (when stable) of the item to the total.
4. Remove the item from the weighing pan, then add the next item and continue as above.
5. Press **RESULT** to view the results from the totalization. (See Figure 63)

Figure 63

Totalization	
Item	Result
Samples	3
Total	506.95 g
Average	168.98 g
$\sigma$ (stdev)	38.90 g
Minimum	117.00 g
Maximum	210.57 g
Range	93.57 g
<div style="display: flex; justify-content: space-around;"> <span>Save to USB</span> <span>Exit</span> </div>	

6. When finished, press the **CLEAR TOTAL** button to reset the accumulated weight to zero.
7. Press **SAVE TO USB** to save the result to a USB flash drive or **EXIT** to return to the totalization home screen.

**NOTE:** Changing units converts the accumulation results to the selected unit.

## FORMULATION

Use this application for compounding and recipe making. (See Figure 64) The maximum number of components is 50.


Figure 64

1. In the upper left portion of the home screen, select **Formulation**. (See Figure 65)
2. Press **START** to begin the compounding process.
3. Place the first ingredient on the pan (or in a tared container) and press **ACCEPT** to store the component.
4. Continue adding components and pressing **ACCEPT** to store the weight of the individual components until the formula is complete. The Total line shows the total weight of all the components.
5. Press **STOP** to finish the formulation. The formulation results are displayed:

**NOTE:** If filler is set to "on" (see Item Settings on page 17), the balance will ask to add a filler material to complete the formulation. Add the filler material and press **ACCEPT** to complete the formulation and display the results.

Figure 65

Formulation	
Comp.	Weight
1	241.76 g
2	272.95 g
Comp. Total	514.71 g
<div style="display: flex; justify-content: space-around;"> <span>Save to USB</span> <span>Exit</span> </div>	

6. Press **SAVE TO USB** to save the results to a USB flash drive or **EXIT** to return to the Formulation Home screen.

**NOTE:** The formulation results will be cleared when a new formulation starts.

## ITEM SETTINGS


To view or adjust the current settings, touch the **ITEM SETTINGS** button. (See Figure 66) The settings screen appears.

**Filler:** If set to "on", a filler material is asked for at the end of the formulation.

**Automatic Mode:** If set to "on", the balance will automatically tare after the component weight has been accepted.

**Print Settings:** Change printing settings.

Figure 66

## WEIGH BELOW

The scale is equipped with a weigh below hook for weighing below the balance.


 **NOTE:** Before turning the balance over, remove the pan and draft shield elements (if present) to prevent damage.

 **CAUTION!** Do not place the balance on the pan support cone or load cell pins.

To use this feature, remove power from the balance, then remove the protective cover for the weigh below opening. (See Figure 67)

The balance can be supported using lab jacks or any other convenient method. Ensure the balance is level and secure. Power on the balance, then use a string or wire to attach items to be weighed.

Figure 67

## MAINTENANCE

### CALIBRATION

Periodically verify calibration by placing an accurate weight on the balance and viewing the result.

### CLEANING

 **WARNING!** Disconnect the balance from the power supply before cleaning.

 **CAUTION!** Do not use solvents, harsh chemicals, ammonia or abrasive cleaning agents.

Make sure that no liquid enters the interior of the balance.

Clean the balance at regular intervals.




Housing surfaces may be cleaned with a lint-free cloth slightly dampened with water or a mild cleaning agent.

Glass surfaces may be cleaned with a commercial glass cleaner. Follow the steps on page 18 on how to remove and install the sliding doors.

## MAINTENANCE CONTINUED

### REMOVING AND REINSTALLING THE GLASS DOORS FOR CLEANING:

1. On the back of the balance, press the pin and slide out the door. (See Figure 68)
2. After cleaning, slide the doors into the slot while pressing the pin mentioned in step one. (See Figure 69)
3. Slide the doors into the slot until the back stopper aligns with the other door. (See Figure 70)


## TROUBLESHOOTING

OPERATING ISSUE	CAUSES	RECOMMENDATIONS
Balance will not turn on.	No power to balance.	Verify connection and voltage.
Poor accuracy.	Improper calibration. Unstable environment.	Platform calibration. Move balance to suitable location.
Cannot calibrate.	Calibration menu locked. Approved Mode set to "on". Unstable environment. Incorrect calibration masses.	Turn Calibration menu lock off. Turn Approved Mode off. Move balance to suitable location. Use correct calibration masses.
Cannot change menu settings.	Sub-menu locked. Approved Mode set to "on".	Unlock sub-menu. Turn Approved Mode off.
Low Reference weight.	Reference weight too small. The weight on the pan is too small to define a valid reference weight.	Increase sample size.
Invalid Piece Weight.	Average piece weight is too small.	Increase average piece weight.
Operation timeout.	Weight reading is not stable.	Move balance to suitable location.
Screen displays "-----".	Busy (tare, zero, printing, waiting for a stable weight).	Wait until completion.

If the troubleshooting section does not resolve your problem, contact Uline Customer Service at 1-800-295-5510.

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
Diagrama 36

Density Determination			
 Density Type	 Auxiliary Liquid	 Water Temp.	 Porous Material
 Print Settings			
Back		Exit	

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Diagrama 37

Density Determination			
 Density Type	 Auxiliary Liquid	 Porous Material	 Print Settings
Back		Exit	

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Density Determination		Max 4200 g d=0.01 g	
<div style="border: 1px solid black; padding: 5px; display: inline-block;"> <p style="font-size: 2em; margin: 0;">1.4541</p> <p style="font-size: 0.8em; margin: 0;">g/cm3</p> </div>			
Dry Weight:	249.94 g	Liquid Weight:	78.37 g
Oil Weight:	249.94 g	Water Temp.:	20.0° C
Item Settings	Start	Menu	

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Diagrama 40

Density Determination			
 Density Type	 Sinker Volume	 Liquid Temp.	 Print Settings
Back		Exit	

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

Diagrama 42

Density Determination		Max 4200 g	d=0.01 g
Weigh sinker in air. Press Accept.			
* <b>210.94</b> g			
Air Weight:	-----	Sinker Volume:	10.0 ml
Liquid Weight:	-----		
Item Settings	Stop	Accept	Menu

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Diagrama 43

Density Determination		Max 4200 g	d=0.01 g
Weigh sinker in liquid. Press Accept.			
* <b>249.94</b> g			
Air Weight:	421.44 g	Sinker Volume:	10.0 ml
Liquid Weight:	-----		
Item Settings	Stop	Accept	Menu

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Diagrama 45

Density Determination			
 Density Type	 Auxiliary Liquid	 Water Temp.	 Porous Material
 Oil Density	 Print Settings		
Back		Exit	

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