

ULINE H-8862

INFRARED FOREHEAD THERMOMETER

1-800-295-5510
uline.com

OVERVIEW OF CONTROLS

CONTROL PANEL PARTS


#	DESCRIPTION
1	Probe
2	Start Button
3	Forehead Mode/Object Mode Button  / 
4	Power On/Off Button 
5	Battery Cover

SAFETY

 **NOTE:** Read all instructions thoroughly before using this product.

The infrared forehead thermometer measures core body temperature, which is the temperature of a body's vital organs. (See Figure 1)

It is specifically designed for safe use on the forehead. It measures body temperature by detecting the intensity of infrared light emitted from the forehead and converts the measured heat into a temperature reading displayed on the LCD. This product is intended for the intermittent measurement of human body temperature by people of all ages. When properly used, it will quickly assess temperature in a timely manner.


WARNING! This product is not intended as a substitute for consultation with a physician.

- Do not allow small children to take their own temperatures unsupervised. There are small parts that can be swallowed.
- Never immerse in water or other liquids.
- Do not modify this equipment without authorization from manufacturer.
- Do not expose thermometer to temperature extremes (below -13°F or over 131°F) or excessive humidity (>95%RH).
- Keep the battery away from children.
- Remove battery from device when not in operation for a long time.

OVERVIEW OF CONTROLS CONTINUED

LCD DISPLAY

FUNCTION	ACTION
Real Time Clock	The real time clock will be recorded with the memory function and help to recognize each measurement result. See the Real Time Clock Setting section on page 3 to learn how to set time functions before first use.
Forehead Mode	Compare the measurement result to regular body temperature. See the Forehead Mode section on page 4 to learn how to measure forehead temperature.
Object Mode	The object mode shows the actual, unadjusted surface temperatures, which is different from the body temperature. It can help monitor the object temperature suitable for the patient. See the Object Mode section on page 4 to learn how to measure object temperature.
Beep Alarm	If thermometer detects a body temperature $\geq 100.0^{\circ}\text{F}$ under forehead mode, there will be a long beep sound followed by three short beep sounds to warn the user.
Memory Mode	Stores 10 readings for forehead and object measurements. Each reading also records the measurement date/time/mode icon.
°C / °F Switch	See Selecting a Temperature Scale section on page 4 to change between celsius and fahrenheit.

REAL TIME CLOCK SETTING

When using thermometer for the first time, set time functions. With thermometer off, press and hold the power on/off button Φ for three seconds to enter setting mode.

SET THE TIME FORMAT

The device can display the time in either an AM/PM (12-hour) or 24:00 (24-hour) format. With the preferred time format on the display, press the the power on/off button Φ . The hour format flashes automatically. (See Figure 2)


Figure 2

SET THE HOUR

Press and release the start button to advance one hour until the correct hour appears. (See Figure 3) After the hour is set, press the power on/off button Φ , the minute figure flashes automatically.

Figure 3

SET THE MINUTE

Press and release the start button to advance one minute until the correct minute appears. After the minute is set, press the power on/off button Φ , the year figure flashes automatically. (See Figure 4)

Figure 4

SET THE YEAR

Press and release the start button to advance one year until the correct year appears. (See Figure 5) After the year is set, press the power on/off button Φ , the month figure appears.


Figure 5

SET THE MONTH

Press and release the start button to advance one month until the correct month appears. After the month is set, press the power on/off button Φ , the date figure flashes automatically. (See Figure 6)

Figure 6

SET THE DATE

Press and release the start button to advance one day until the correct day appears. (See Figure 7) After the day is set, press the power on/off button Φ to exit the setting mode.


Figure 7

TEMPERATURE-TAKING HINTS

To ensure that the reading always reflects body temperature accurately, note the following factors that may affect an accurate reading:

1. Know each individual's normal temperature when they are well. This is the only way to accurately diagnose a fever. To determine normal temperature, take multiple readings when healthy. Re-measure with a standard digital thermometer for confirmation.
2. Users and the thermometer should be in the same ambient temperature for at least 10 minutes before taking a reading.
3. Users should not drink, eat or be physically active such as bathing, showering, shampooing and hair drying before/while taking the measurement.
4. Oils or cosmetics on the forehead may give a lower temperature reading than the actual one. Remove dirt from the forehead before taking a measurement.
5. Do not hold a hand on the forehead for any length of time.
6. Do not take temperature over scar tissue, open sores or abrasions.
7. Do not use the thermometer on a perspiring or sweating forehead.
8. Do not take a measurement while or immediately after nursing a baby.
9. Thermometer is intended for home and clinical use.
10. Do not take temperatures near places that are very hot, such as fireplaces and stoves.
11. Do not touch the probe window. The probe window of the thermometer is the most delicate part of the device. The accuracy of the reading may be affected if the probe window is damaged or dirty.
12. If thermometer is stored in a significantly different environment than testing location, place it in testing location for approximately 30 minutes prior to use.
13. It is not intended for use in oxygen rich environments or in the presence of flammable anesthetic mixture with air, oxygen or nitrous oxide.

SELECTING A TEMPERATURE SCALE

1. Temperature readings are available in the Celsius (°C) or Fahrenheit (°F) scale. (See Figure 8)
2. With the thermometer off, press and hold the start button for three seconds to enter unit changing mode.
3. Press and release the start button to select the unit.


Figure 8

When the preferred unit appears on the display, press the power on/off button ϕ to exit unit changing mode.

CHANGING MODES

Press the forehead mode/object mode button $\text{person}/\text{house}$ to switch between forehead mode and object mode.

FOREHEAD MODE

1. Press the power on/off button ϕ . The display is activated to show all segments. After self-checking, the start up screen appears on the display screen with beeps, indicating that a new measurement can be started. (See Figure 9)
2. Aim the thermometer at the center of the forehead with a distance of 2-3 cm and then press the start button. (See Figure 10)

Figure 9

Figure 10

 **NOTE:** Do not remove the thermometer from the forehead before hearing the final beep. A waiting period of three seconds between testing is recommended to ensure complete accuracy.

3. Read the temperature on the display.
4. Press the power on/off button ϕ to turn off.

BACKLIGHT

1. The display will turn green for three seconds with a happy face ☺ when the unit is ready for measurement and when a measurement is completed with a reading less than 99.1°F.
2. The display will turn yellow for three seconds with a happy face ☺ when a measurement is completed with a reading less than 100°F.
3. The display will turn red for three seconds with a sad face ☹ when a measurement is completed with a reading equal to or higher than 100°F.

OBJECT MODE

1. Press the power on/off button ϕ to turn on the thermometer. User can take the object temperature after two beep sounds. (See Figure 11)
2. Aim the thermometer at the center of the object to measure with a distance of 1-2 cm.
3. Press the start button, and then read the temperature on the display.
4. Press the power on/off button ϕ to turn off.


Figure 11

AFTER MEASUREMENT

1. Device will automatically shut off if left idle for more than one minute to extend battery life.
2. Clean probe after each use to ensure an accurate reading and prevent cross-contamination. See the Care and Cleaning section on page 5 for details.

MEMORY MODE

The memory mode can be accessed in either mode.

1. When the thermometer has been turned on or finished testing, press and hold the power on/off button ϕ for three seconds. The letter M will appear in the center right corner of display. (See Figure 12)
2. The thermometer will automatically memorize the last 10 temperature readings. Each memory also records the measurement date/time/mode icon. Each time the power on/off button ϕ is pressed, the screen displays past readings that correspond with a number 1-10. Number 1 reflects the most recent reading, while number 10 reflects the oldest reading stored in memory. (See Figure 13)
3. In memory mode, the forehead mode mark person and object mode mark house always exist. User can press the start button to take new measurements.


Figure 12

Figure 13

CARE AND CLEANING

1. Probe window must be kept clean, dry and undamaged at all times. Accuracy of temperature readings can be affected by damage to the probe window, the presence of dirt, fingerprints, earwax, dust and other soiling compounds on the probe window. Degraded sensors can reduce performance or cause other problems.
2. Periodic cleaning and disinfection of the device can prevent patient cross-infection. Clean and disinfect thermometer immediately after use.

CLEANING

1. Soak a clean, soft cloth in drinking water, wring it out dry and then wipe the thermometer (including probe) no fewer than three times.
2. Use another clean soft cloth to wipe the thermometer dry.
3. Put the thermometer in the original packaging.

DISINFECTING

1. Soak a clean, soft cloth in drinking water, wring it out dry and then wipe the thermometer (including probe) no fewer than three times.
2. Use another clean, soft cloth to wipe the thermometer dry.

3. Using a clean, soft cloth dipped in 70% medical alcohol, wipe the probe three times, at least one minute for each time.
4. Using a clean cotton swab dipped in 70% medical alcohol, wipe the sensor window three times.
5. Wait at least 10 minutes to let the alcohol evaporate and then put the thermometer in the original packaging.
6. Use a soft, dry cloth to clean the thermometer display and exterior.
7. The thermometer is not waterproof. Do not submerge the unit in water when cleaning.
8. Store in a dry location, free from dust and contamination and away from direct sunlight.
9. Strong electromagnetic fields may interfere with the proper operation of thermometer. The device needs special precautions regarding EMC according to the EMC information.
10. Store thermometer in original packaging after use.

BATTERY REPLACEMENT

1. Replace battery when the battery low symbol  appears in the upper right corner of LCD display. (See Figure 13)
2. Slide battery cover down as shown. (See Figure 14)
3. Remove battery and install two new AAA alkaline batteries.
4. Slide battery cover back on. (See Figure 15)

Figure 13


Figure 14


Figure 15

SPECIFICATIONS

REFERENCE	DESCRIPTION
Measuring Range	Forehead Mode: 93.2°~109.4°F (34.0°~43.0°C) Object Mode: 32°~212°F (0°~100°C)
Measuring Site	Forehead (Forehead Mode)
Reference Body Site	Oral (This thermometer converts the forehead temperature to display its "oral equivalent".)
Operation Mode	Forehead Mode (Adjust Mode)
Laboratory Accuracy	Forehead mode: ± 0.4°F (0.2°C) during 95.9°~107.6°F (35.5°~42.0°C) at 59.0°~95.0°F (15°~35°C) operating temperature range ± 0.5°F (0.3°C) for other measuring and operating temperature range Temperature range: Object mode: ± 4% or ± 4°F (2°C) whichever is greater
Display Resolution	0.1°F or 0.1° C
Measure Time	Approximately one second
Operating Temperature Range	50°~104°F (10°~40°C), 15%~85%RH, non-condensing Atmospheric Pressure: 700~1,060hPa
Storage and Transport Temperature Range	-13°~131°F (-25°~55°C), 15%~85%RH, non-condensing Atmospheric Pressure: 700~1,060hPa
Clinical Accuracy	Clinical Bias: -0.16°F (-0.09°C) Clinical Repeatability: 0.23°F (0.13°C) Limits of Agreement: 1.57°F (0.87°C)
Shock	Withstands drop of three feet
Dimension	153 x 41 x 44 mm
Weight	Approximately 84 grams (with batteries)
Battery	DC3V (2XAAA battery)
Battery Life	Approximately 1 year/6,000 readings
Expected Service Life	Three years
Ingress Protecting Rating	IP22

TROUBLESHOOTING



ERROR MESSAGE	PROBLEM	SOLUTION
Er 1	Measurement before thermometer is ready	Take a measurement until the forehead mode icon  or object mode icon  appears on the display.
Er 2	The ambient temperature is not within the range between 50°~104°F.	Place thermometer in a room for at least 30 minutes at room temperature between 50°~104°F.

TROUBLESHOOTING CONTINUED

ERROR MESSAGE	PROBLEM	SOLUTION
Er3	The thermometer is placed incorrectly or unsteady.	Read instructions thoroughly and take a new temperature measurement.
Er4	The thermometer is showing a rapid ambient temperature change.	Allow the thermometer to rest in a room for at least 30 minutes at room temperature: between 50°~104°F.
Er5	The thermometer is not functioning properly.	Unload the battery, wait for one minute and repower it. If message appears, contact Uline Customer Service at 1-800-295-5510.
Hi	<p>In Forehead Mode: Temperature taken is lower than 109.4°F.</p> <p>In Object Mode: Temperature taken is lower than 212°F.</p>	Read Temperature Taking Hints thoroughly (see page 3), then take the new temperature measurement.
Lo	<p>In Forehead Mode: Temperature taken is lower than 93.2°F.</p> <p>In Object Mode: Temperature taken is lower than 32°F.</p>	Read Temperature Taking Hints thoroughly (see page 3), and ensure the lens and filter are clean. Take the new temperature measurement.
	The thermometer does not work due to low battery.	Replace old batteries with two new AAA alkaline batteries.

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