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
1-800-295-5510
uline.com

SAFETY

Balance must only be used for measuring the parameters described in these instructions. Any other type of use and operation beyond the limits of technical specifications, without written consent from the manufacturer, is prohibited.

This instrument complies with current industry standards and the recognized safety regulations; however, it can constitute a hazard in use. If the balance is not used according to these instructions, the intended protection may be impaired.

CAUTION! Read all safety warnings before installing, making connections or servicing the balance. Failure to comply with these warnings could result in personal injury and/or property damage. Keep instructions for future reference.

- Before connecting power, verify that the power adapter's input voltage range and plug type are compatible with the local AC mains power supply.
- Position the balance in a location where power connection is easily accessible. Ensure the power cord does not pose a potential obstacle or tripping hazard.
- Operate the balance only under ambient conditions specified in these instructions.
- Balance is for indoor use only.
- Do not operate the balance in wet, hazardous or unstable environments.
- Do not allow liquids to enter the balance.
- Do not load the balance above its rated capacity.
- Do not drop loads on the platform.
- Do not place the balance upside down on its platform.
- Use only approved accessories and peripherals.
- Disconnect the balance from the power supply when cleaning.
- Service should only be performed by authorized personnel.
- When shipping or transporting the balance, follow the applicable regulations for equipment containing lithium-ion batteries.

WARNING! Never work in an environment subject to explosion hazards. This instrument's housing is not gas tight. The ingress of gases may cause corrosion and spark formation, resulting in a risk of explosion.

WARNING! Electrical shock hazards! The housing should only be opened by authorized and qualified personnel. Remove all power connections to the unit before opening.

CAUTION! The protection provided by the unit may be impaired if used with accessories not provided or recommended by the manufacturer or used in a manner not specified by the manufacturer.

LOCATION

External factors such as temperature fluctuations, air currents, electromagnetic interference and vibrations can all have an impact on accuracy.

Keep the balance in a controlled environment to minimize these factors and ensure accurate readings.

INSTALLATION

Unpack the balance and install the weighing pan. Save the packaging for safe storage and transportation.

POWER

1. Connect the DC output connector to the power receptacle on the rear of the base, then connect the plug to a suitable electrical outlet. (See Figure 1)

CAUTION! For use with CSA-certified (or equivalent approved) power source, which must have a limited current output.

CAUTION! Allow equipment to warm up for 60 minutes for optimal weighing performance.


Figure 1

2. Press the ***POWER ON*** button on the front of the balance to turn it on. (See Figure 2)


Figure 2 Power On Button

3. Log in to the balance with default user ID and password. (See Figures 3-6)

NOTE: User ID: Admin; Password: No password required

Figure 3

Figure 4

Figure 5

Figure 6


4. When the balance is in sleep or standby mode, press the ***CLICK HERE*** icon to start weighing. (See Figure 3)

INITIAL LEVELING

Balance is equipped with a digital leveling bubble. The digital level indicator in the upper left corner will notify user if the weighing table is not level.

1. Press the leveling button. (See Figure 7)

Figure 7

2. Adjust the leveling feet according to the position of the digital leveling bubble until the bubble is centered. (See Figures 8-9)

NOTE: If the environment is unstable or the weighing table is tilting too much, the leveling process may fail. In such cases, move the balance to a stable, level surface that can absorb vibrations.

Figure 8



Figure 9


REMOVING DRAFT SHIELD

NOTE: For H-12947 (EXR224) only.


1. Remove the weighing pan. Turn the fixing ring clockwise to unlock the EMC plate. (See Figure 10)

Figure 10


2. Remove the EMC plate and the four screws beneath it. Lift the draft shield vertically to remove it from the scale body. (See Figure 11)

Figure 11


H-12947 (EXR224) FEATURES

FRONT VIEW


BACK VIEW


H-12947 (EXR224) FEATURES CONTINUED

BOTTOM VIEW


H-12948 (EXR4202) AND H-12949 (EXR6202) FEATURES

FRONT VIEW


BACK VIEW

OVERVIEW OF CONTROLS

CONTROL PANEL

#	DESCRIPTION
1	Upper Left Shortcut Buttons
2	Top Line
3	Upper Right Icons
4	Menu Button
5	Main Display Area
6	Lower Left Icon
7	Touchless IR Sensors
8	Power On Button
9	Item Settings Button
10	Zero Button
11	Print Button
12	Tare Button
13	Calibration Button
14	Library Button
15	Add+

UPPER LEFT SHORTCUT BUTTONS

- Leveling
- Batch Print

TOP LINE

Displays essential balance information, such as Model Name, Max Capacity, Displayed Readability Value (d), Verification Scale Value (e), User Name, Date/Time.

UPPER RIGHT ICONS

- Peripheral icons connected to balance.
- Wi-Fi Connectivity
- Real-time Leveling Bubble
- Warning Indicator

MENU BUTTON

Switch between main menus.

MAIN DISPLAY AREA

- Click "g" to switch between weighing units.
- Reference field shows parameters for the current weighing application.

LOWER LEFT ICON

Change weighing application mode.

HOT KEYS

- Item Settings, Zero, Print, Tare, Calibration, Library
- Press ***ADD***+ to add customized hot keys.

ADD+

- 1d/10d
- Digital Leveling
- Repeatability Test

HOT KEYS

Bottom of main screen has eight hot keys for basic weighing applications. Customize an additional five hot keys by using the ***ADD+*** button.

BUTTON	FUNCTION
ITEM SETTINGS	Configure current application mode settings.
ZERO	Remove the load from the pan and press ***ZERO*** to set the display to zero. When weighing pan is empty, the >0< indicator turns on when the measurement is within +/– 1/4 division (d) of the zero setting. **NOTE: The balance features Auto Zero Tracking (AZT) that automatically maintains a center of zero within the tolerances set in the AZT menu. See Auto Zero Tracking on page 12.**
PRINT	Send the displayed value to a printer or computer via the active COM port. Ensure the balance has connected with printer or PC and interface parameters are set up correctly.
TARE	Tare the weight of the container to zero, and the net weight will show on the reference field.
CAL	Calibrate the balance or set the calibration parameters.
LIBRARY	Create a new library or retrieve an existing one for the current application.
SENSORS	The balances have two pressless sensors that can be assigned a unique function when activated. The Sensor setting are: Off, Print, Zero, Calibration, and Tare. **NOTE: Once the sensor is set up for a specific function, waving your hand over it will trigger that function and cause the sensor light to turn green. If the sensor fails to activate, verify the settings. The sensor lights will remain red.**
SAVE TO USB	When the user inserts the USB flash drive, an icon will appear on the right side of the Sensor. While data is being saved to the USB drive, the icon will display the progress. **NOTE: Do not unplug the USB drive during data transfer.**
ADD+	Access five customizable hot keys. Click "Right Arrow" key to move the available hotkey to selected area and press "Save" to exit. The new hotkey will appear on the right side of the "Sensor" section. Swipe along the bottom line, you will find the new hotkey.

ADDITIONAL HOT KEYS

1. Click right arrow key to move the available hot key to selected area. Press ***SAVE*** to exit.
2. The new hot key will appear on the right side of the ***SENSOR*** section. Swipe along the bottom line to find the new hot key.

BUTTON	FUNCTION
1d/10d	Change division for weighing result.
Digital Leveling	Level the balance.
Repeatability Test	Perform the repeatability test.

BALANCE MENUS

MENU NAVIGATION

1. **Enter the Menu:** The main menu button is hidden on the middle-left side. Press the button to reveal the menu.
2. **Select the Sub-Menu:** Locate the item on the main menu list and touch it. The display highlights the item in red for about one second before the sub-menu appears.

MAIN MENU	SUB-MENUS
Quick Setup	Language, Time Synchronization, Time Zone (UTC), Date and Time, Automatic Calibration, User Management, System Log, Balance Info, Digital Leveling, Repeatability Test
Calibration	Internal Calibration, Automatic Calibration, Span Calibration, Perform Initial Calibration, Perform Span External Calibration, Calibration History
Balance Setup	Language, Time Synchronization, Time Zone (UTC), Date and Time, Balance Name, Change Password, Fingerprint, Fingerprint Setting, Stability Indicator Range, Filter Level, Auto Zero Tracking, Gross Indicator, Graduation, Ionizer, Auto Doors, Sensors, System Log, ECO
User Management	User Management, Group User Permissions, Password Policy
Application Modes	Basic Weighing, Parts Counting, Check Counting, Percent Weighing, Check Weighing, Dynamic Weighing, Totalization, Formulation, Differential, Peak Hold, Density Determination
Weighing Units	Vary by model. See Model Specifications on page 55.
Communication	RS-232, USB, Ethernet, Wi-Fi and Bluetooth
Library	Library List, Delete Library, Import Library, Export Library
Maintenance	Diagnostic Menu, Software Upgrade, Service Log, Service Menu
Factory Reset	Reset All, Quick Setup, Calibration, Balance Setup, Application Modes, Weighing Units, Communication, Library

QUICK SETUP

The quick setup menu is designed for first-time users. Navigate through this menu to configure desired functions. For details on each setting, see Balance Setup on page 12.

- **Language:** Set the language for menu display and pop-up messages. The default setting is configured based on the country of sale.
- **Time Synchronization:** Balance features synchronization with network (NTP system), which can use consistent time and data from the local network. The NTP system supports DNS domain name resolution.
- **Date and Time:** Set UTC time zone, time format and date format.
- **Automatic Calibration:** Enable or disable AutoCal™. When enabled, AutoCal™ will automatically calibrate the balance using the internal mass.
- **User Management:** Create and manage up to 200 user accounts, featuring four predefined roles. See User Management on page 13.
- **System Log:** Balance can store up to 100,000 log entries in four log options. See System Log on page 13.
- **Balance Info:** Display will show information about the balance.
- **Digital Leveling:** The digital level bubble will indicate level status of the balance.
- **Repeatability Test:** This function conducts the daily repeatability test of the balance. Users can perform a ten-time repeatability test using weights that are close to their daily working point.

REPEATABILITY TEST PROCEDURE

NOTE: To determine the minimum weight in accordance with USP41, test environment must be equivalent to standard lab USP1205, and the test weights must be certified and their values checked periodically.

1. Press ***REPEATABILITY TEST*** and follow the on-screen instructions.
2. In repeatability test settings, set Number of Tests (default is 10), Test Weight ID, Test Weight Value and Test Weight Class in accordance with the test environment's standard procedure.
3. Press ***SAVE*** to store the settings.
4. Place the balance in a stable environment and press ***ZERO*** or ***TARE*** before starting.
5. Press ***START*** and place weights on the pan one by one. During the test, do not tare or zero the pan, as the load will calculate all factors to produce accurate results.
6. After set number of tests has been completed, the balance will automatically display the results.
7. Press ***VIEW HISTORY*** to review the last repeatability test result.
8. Press ***PRINT*** to print the repeatability test report and export it to a PDF file when a USB flash drive is inserted.
9. Press ***MIN-WEIGHT*** to display the reference minimum weight.

CALIBRATION

The balance offers three calibration methods: Automatic Calibration, Internal Calibration and Span Calibration. Calibration is accomplished with the internal calibration mass and can be performed at any time, provided the balance has warmed up to operating temperature and is level.

NOTE: Do not disturb the balance during calibration.

CALIBRATION SETTINGS

- **Internal Calibration:** On or off. Default setting is on.
- **Automatic Calibration:** On or off. Default setting is on.
- **AutoCal™ Interval (Hours):** Set interval between automatic calibrations. Settings of 3-11 hours are allowed.
- **Print Automatic Calibration:** On or off. Default setting is on.
- **Span Calibration:** Set span calibration points. See Span Calibration on page 11.

INTERNAL CALIBRATION

1. Press ***INTERNAL CALIBRATION***, then press ***START***. The balance will begin to calibrate. Press ***CANCEL*** at any time to stop calibration.
2. After calibration, place the test mass on the pan and verify that the mass value matches the displayed value. If not, repeat steps until the reading agrees with the test mass.

CALIBRATION REPORT

- When the calibration process is completed, the balance will generate a calibration report that includes the following information:
 - **Date/Time**: Week, Date and Time.
 - **Balance ID:** Balance series number.
 - **Balance Name:** Balance Model.
 - **Result:** Successful or Failed.
 - **Signature:** User ID (when the user management is set on).
 - **Verified by:** User can input the User ID/ Password.
- Press the Print button to print the calibration report, and the results will be saved in the Calibration History for future reference.

AUTOMATIC CALIBRATION

When internal calibration is enabled, balance will perform a self-calibration:

- When environment temperature changes by 34.7°F.
- Every 11 hours.
- An information screen will appear when an automatic calibration is about to start. Three option buttons will be displayed:
- **Now:** Press to perform the calibration immediately.
- **After 5 min:** Press to perform the calibration after five minutes.
- **Deactivate:** Press to deactivate the automatic calibration function.

AutoCal™ Interval (Hours)

Automatic calibration will be triggered on a time-based schedule.

1. Enable automatic calibration.
2. Set number of hours between each automatic calibration. Value can be set from 3 to 11 hours.

Print Automatic Calibration

When enabled, the balance will send the calibration report to a connected PC or printer. Results will be saved in calibration history for future reference. The report will include the following information:

- **Date/Time:** Week, Date and Time
- **Balance ID:** Balance Series Number
- **Balance Name:** Balance Model
- **Result:** Successful or Failed
- **Signature:** User ID (when user management is enabled)
- **Verified by:** User can input the User ID/Password.

SPAN CALIBRATION

User can choose the base span calibration points by using the full load or half load of the balance's capacity. Define span calibration weights by preset values or alternative weights (between 20 and 100% of the maximum balance capacity).

MODEL	SPAN CALIBRATION POINTS	WEIGHT CLASS	
H-12947 (EXR224)	50 g, 100 g, 150 g, 200 g*	ASTM Class 1	OIML E2
H-12948 (EXR4202)	1,000 g, 2,000 g, 3,000 g, 4,000 g*	ASTM Class 1	OIML E2
H-12949 (EXR6202)	2,000 g, 3,000 g, 4,000 g, 6,000 g*	ASTM Class 1	OIML E2

*Denotes default setting

1. Press ***SPAN CALIBRATION***, then press ***START***. The balance will begin to calibrate. Press ***CANCEL*** at any time to stop calibration.
2. Follow the guide message to place weights on the pan and wait until the balance is stable.
3. After calibration, place the test mass on the pan and verify that the mass value matches the displayed value. If not, repeat steps until the reading agrees with the test mass.

Print Span Calibration

- When enabled, the balance will send the calibration report to a connected PC or printer. Results will be saved in calibration history for future reference. The report will include the following information:
- **Date/Time:** Week, Date and Time
- **Balance ID:** Balance Series Number
- **Balance Name:** Balance Model
- **Result:** Successful or failed, and difference deviation compared with the latest calibration value
- **Signature:** User ID (when user management is enabled)
- **Verified by:** User can input the User ID/Password.

CALIBRATION HISTORY

Balance can store up to 3,000 calibration histories and log when calibration is performed. The total memory is displayed as xx/xx% in the upper right corner. Calibration records can be exported to a PDF file when a USB flash drive is inserted.

NOTE: When user activates the system log function, calibration history will be recorded as a calibration log in the service menu. See System Log on page 13.

BALANCE SETUP

Customize balance functionality.

LANGUAGE

Set the language for menu display and pop-up messages. The default setting is configured based on the country of sale.

Available Languages: English, German, French, Spanish, Italian, Polish, Czech, Hungarian, Portuguese, Chinese, Japanese, Korean, Turkish

TIME SYNCHRONIZATION

Balance features synchronization with network (NTP system) which can use consistent time and data from the local network. The NTP system supports DNS domain name resolution.

When user inputs the address of a public domain server, the balance will automatically convert it into an IP address. For example, the US would use "0.us.pool.ntp.org".

DATE AND TIME

Set UTC time zone, time format and date format.

- **Time Formats:** 12 Hours, 24 Hours
- **Date Formats:** YYYY-MM-DD, MM-DD-YYYY, DD-MM-YYYY

BALANCE NAME

Set the balance ID with up to 25 alphanumeric characters. Default setting is the model name of the balance.

CHANGE PASSWORD

Change the password for the user currently logged in. See User Management on page 13.

- **Alphanumeric:** Password can contain 8-10 characters, combining letters and numbers.
- **Numeric:** Password can contain up to 25 characters.

NOTE: If user forgets their password, contact Uline Customer Service at 1-800-295-5510.

FILTER LEVEL

Set the amount of signal filtering. Higher numbers mean slower stabilization time with more stability.

Available Settings: 1, 2, 3, 4, 5. Default setting is 5.

STABILITY INDICATOR RANGE

Set the amount the reading can vary while the stability symbol remains on. Number corresponds to stability factor.

Available Settings: 0.25, 0.5, 1, 2, 3, 4. Default setting is 1.

AUTO ZERO TRACKING

Set the automatic zero tracking functionality. Example: 0.5d = display maintains zero up to a drift of 0.5 graduation per second.

Available Settings: Off, 0.5d, 1d, 3d

GROSS INDICATOR

Set the symbol displayed for gross weights.

Available Settings: Off, gross (G), brutto (B)

GRADUATION

Set the displayed readability of the balance. Example: 0.01 g readability means selecting 10 divisions will result in a displayed reading of 0.1 g.

Available Settings: 1d (standard), 10d (readability increased by a factor of 10)


SENSORS

Two touchless IR sensors are located on top of the balance terminal. When user waves a hand over the sensor, the function will be activated and a green light will flash. (See Figure 12)

Available Settings: Off, Print, Zero, Tare, Calibrate

NOTE: If sensor is not activated, the sensor light will be red.

Figure 12

SYSTEM LOG

All changes made to the balance settings will be recorded in the system log file (printed data, date/time changes, setting changes, calibration actions, user log ins/log outs, user account creation/edits/deletion, etc.).

Internal system log capacity is 100,000 entries. When memory is full, pop-up will notify user to export log files to a USB flash drive.

System log is saved in the service menu.

- **User Management:** Create and manage up to 200 user accounts, featuring four predefined roles. See User Management below right.
- **Auto Login:** Beneficial for admin users who have not set a password. No login is required after pressing the ***STANDBY*** button. Other users cannot use this function when user management is enabled.
- **Electronic Signature:** When enabled, User ID will be used in the printout field.
- **System Log:** Enable or disable. When disabled, changes to the balance settings will not be recorded in balance memory.

ECO

Set energy- and power-saving settings.

Power Saving

- **Off:** Touch screen display will operate with default settings.
- **Auto Dim:** Display will shift to a darker shade when inactive.
- **Auto Standby:** Balance will automatically log off.
- **Sleep Mode:** Display will be turned off at a predetermined time, only retaining base operation power.

ECO Time

- **Auto Dim or Auto Standby:** ECO time will automatically set to 20 minutes.
- **Sleep Mode:** ECO time will automatically set to 0 minutes.
- For all settings, brightness, volume and terminal status light will be locked to conserve power.

Brightness

Adjust screen brightness.

Available Settings: 0-100. Default setting is 90.

Volume

Balance has sound options for clicks, successful actions and errors.

Available Settings: 0-20. Default setting is 2.

Terminal Status Light

Customize lighting system.

Available Settings: Off, Red, Pink, Yellow, Green, Cyan, Blue, White

USER MANAGEMENT

When enabled, the balance will be able to create up to 200 user accounts, featuring four predefined roles.

- **Administrator:** 1
- **Supervisor:** Maximum of 20
- **Operator, Log Viewer, Group 1 Users**: 179

NOTE: User ID is generated from first and last name entries. User ID is unique and cannot be changed once it is confirmed.

CREATE, EDIT OR DELETE A USER

- **New**: Add a new user account
- **Edit**: Modify an existing user account.
- **Delete and Delete All**: Remove user accounts. Only available to Administrators and Supervisors.

NOTE: Administrator cannot be deleted.

PREDEFINED USER ROLE ACCESSIBILITY

ROLE	ADMINISTRATOR	SUPERVISOR	OPERATOR	LOG VIEWER	GROUP 1 USER
Number of Users	1	20	Total of 179		
Create Account/Delete Account	Create/Edit/Delete Supervisor, Operator, Log Viewer and Group 1 User	Create/Edit Operator	X	X	Access permissions depend on user settings. See Group User Permissions on page 15.
Date and Time	√	X	X	X	
Print Content	√	√	X	X	
Calibration	√	√	X	X	
Weighing Units	√	√	√	X	
Library	√	√	New, Load, Edit/Delete (own)	X	
Factory Reset	√	√	X	X	
Repeatability Test	√	√	Allowed to operate, not allowed to change weight settings	X	
Balance Name	√	X	X	X	
Set Up a Min-Weight	√	√	X	X	
Service Mode	√	X	X	X	
System Log	Read/Export/Delete	Read/Export	X	Read/Export	
Graduation 1d/10d	√	√	X	X	

GROUP USER PERMISSIONS

Allows for the streamlined assignment of identical access permissions to multiple users.

- **Predefined Groups:** Administrator, Supervisor, Operator, Log Viewer
- **Customizable Group 1:** Configure unique access permissions.

All access options, except for Date and Time, Print Content, and Balance Name, are pre-configured and ready for use.

PASSWORD POLICY

Allows password protection for each user.

- **Expiration Days:** 30 days, 60 days, 90 days
- **Password Requirements:** Numeric (up to 25 characters), Alphanumeric (8-10 characters, combining letters and numbers)

NOTE: If user forgets their password, contact Uline Customer Service at 1-800-295-5510.

Press ***LOG OFF*** to log out of the current user account.

APPLICATION MODES

Balance can be configured to operate in 11 application modes: basic weighing, parts counting, check counting, percent weighing, check weighing, dynamic weighing, totalization, formulation, differential, peak hold, density determination.

NOTE: If an application mode does not appear in the list of available applications, it must be enabled in the main menu.

1. Touch ***MENU***, then select "Application Modes". (See Figures 13-14)

Figure 13




Figure 14

2. The full list of modes appears. Enable or disable application modes as desired. Press ***BACK*** to exit to main screen. (See Figure 15)

Figure 15

3. To switch modes, press active weighing application mode in lower left side of display.

NOTE: Before using any weighing application, ensure balance has been leveled and calibrated successfully.

WEIGHING UNITS


Balance can be configured to measure in a variety of units, including two customized units. The default unit is grams.

STANDARD UNITS

On the main screen, click "g" to switch to an alternative weighing unit. (See Figure 16)


NOTE: If a weighing unit does not appear in the list of available applications, it must be enabled in the weighing units menu.

Figure 16

1. Touch ***MENU***, then select "Weighing Units".
2. The full list of weighing units appears. Enable or disable units as desired. (See Figure 17)

Figure 17

3. Press ***BACK*** to exit to main screen.

CUSTOMIZED UNITS

Display weight in an alternative unit of measure. The customized unit is defined using a conversion factor, where the conversion factor is the number of customized units per kilogram expressed in scientific notation (Factor x $10^{Exponent}$).

1. Touch ***MENU***, then select "Weighing Units".
2. Press ***CUSTOMIZED UNITS***.
3. Set unit name, factor, exponent and least significant digit.
4. Turn customized unit on or off.
5. Press ***BACK*** to exit to main screen.

The following example shows how to display weight in troy ounces (32.15075 troy ounces per kilogram):

1. Set Factor of 0.321508.
2. Set Exponent of 2.
3. Set Customized Units status on.

Factor

Set the conversion factor using the numeric keypad.

Available Settings: 0.00001 to 1.9999999. Default setting is 1.000000.

Exponent

Set the factor multiplier.

0 = multiply the factor by 1 (1 x 100)

1 = multiply the factor by 10 (1 x 10^{1})

2 = multiply the factor by 100 (1 x 10^{2})

3 = multiply the factor by 1,000 (1 x 10^{3})

-1 = divide the factor by 10 (1 x 10^{-1})

-2 = divide the factor by 100 (1 x 10^{-2})

-3 = divide the factor by 1,000 (1 x 10^{-3})

Least Significant Digit

Set the graduation.

Available Settings: 0.00001, 0.00002, 0.00005, 0.0001, 0.0002, 0.0005, 0.001, 0.002, 0.005, 0.01, 0.02, 0.05, 0.1, 0.2, 0.5, 1, 2, 5, 10, 20, 50, 100, 200, 500, 1,000


NOTE: Least significant digit available settings are dependent on factor and exponent settings. Selections are limited.

BASIC WEIGHING

Used to determine the weight of items in the selected unit of measure. Basic weighing application mode is displayed by default on the main screen.

1. Select "Basic Weighing" in the lower left side of the display. (See Figure 18)

Figure 18

2. Press ***TARE*** or ***ZERO*** to start weighing.
3. Place sample on the pan. When stable, the * appears.
4. The resulting value is displayed in the main weighing line in the active unit of measure.


REFERENCE FIELD

- **Gross:** Displays the gross weight when the container is tared.
- **Tare:** Displays the tare value when the container is tared.
- **Min-Weight (Conditional):** Displays minimum weight value when the minimum weight is set to "On" or "Below Minimum Weight".

ITEM SETTINGS

Press the ***ITEM SETTINGS*** button to change the application settings. (See Figure 19)


Figure 19

Primary Weighing Unit

Choose weighing unit to display on the main screen. Default setting is grams. (See Figure 20)

Figure 20

Secondary Weighing Unit

When selected, secondary unit will display below the primary unit on the main screen.

Minimum Weight

1. Choose minimum weight setting.
 - **Off:** Disable minimum weight function.
 - **On:** When the result is lower than minimum weight value, the screen will show "Below Minimum Weight", and user will have to add more samples. Result will not appear on printouts.
 - **Below Minimum Weight:** When the result is lower than minimum weight value, the screen will show "Below Minimum Weight", and user will have to add more samples. Result will appear on printouts. (See Figures 21-22)

Figure 21




Figure 22

2. Press ***SAMPLE LIBRARY*** to input the minimum weight value. Set between 0 g and maximum capacity of balance. If number is out of range, "Invalid Entry" message will appear. (See Figure 23)

Figure 23

Auto Tare

Choose whether to automatically tare the container value. Default setting is off.

When auto tare is enabled, the display shows the message "Place container on the pan". After a container is placed on the weighing pan, its weight will be stored in the balance. The net weight and gross weight values will then be displayed in the reference field. (See Figures 24-25)

Figure 24



Figure 25

LIBRARY AND PRINTING SETTINGS

NOTE: For details, see Library on page 45 and Printing Settings on page 46.

PARTS COUNTING

Used to count samples of uniform weight. Displays the number of pieces or parts based on an average piece weight.

NOTE: Set the sample size and average piece weight (APW) before starting parts counting application.

1. Select "Parts Counting" in the lower left side of the display.
2. Press ***TARE*** or ***ZERO*** to start weighing.
3. Place objects on the pan. When stable, the * appears and number of pieces is shown on the display. The default (or last) APW is used.

REFERENCE FIELD

- **Gross:** Displays the gross weight when the container is tared.
- **Net:** Displays the net weight when the container is tared.
- **Tare:** Displays the tare value when the container is tared.
- **APW:** Displays the average piece weight.
- **Sample Size:** Displays the set sample size.
- **Sample Name:** Displays the set sample name.

APPLICATION BUTTONS

BUTTON	FUNCTION
SAMPLE SIZE	Input sample size.
SET APW	Place the sample on the pan. Balance will calculate average piece weight (APW).
INPUT APW	Input an established APW.
0	Updates to display the number of parts counted. Press to review the statistical results.
+	Press to manually accumulate weight. In automatic statistics mode, this button will be inactive.

ITEM SETTINGS

Press the ***ITEM SETTINGS*** button to change the application settings.

Auto Optimization

Automatically optimize APW to improve counting accuracy.

Statistics

Choose statistics mode.

- **Automatic:** Weights are automatically recorded when stable.
- **Manual:** User must press the **+** button to record statistical data.

Secondary Weighing Unit

Choose between pre-set weighing units and two customized units.

Auto Tare

Choose whether to automatically tare the container value. Default setting is off.


When auto tare is enabled, the display shows the message "Place container on the pan". After a container is placed on the weighing pan, its weight will be stored in the balance. The net weight and gross weight values will then be displayed in the reference field.

SET A KNOWN AVERAGE PIECE WEIGHT

Each time a new type of part is counted, the nominal weight of one piece must be established using a small quantity of pieces.


1. Press ***INPUT APW*** and enter the desired APW value. (See Figure 26)

Figure 26

2. Press ***ENTER*** to save the value. (See Figure 27)

Figure 27

SET A DERIVED AVERAGE PIECE WEIGHT

1. Press ***SAMPLE SIZE*** and enter the desired sample size. The sample size can be 1-10,000 pieces. Default is 10.
2. Press ***ENTER*** to save the value.
3. Place reference weight on the pan and press ***ACCEPT*** to establish a new APW. The home screen will display the new APW. (See Figures 28-29)

Figure 28



Figure 29

NOTE: Once a sample size is changed, the balance will immediately recalculate APW screen.

ACCUMULATE DATA

Utilize the Statistics "+" function to aggregate samples based on counting numbers and generate statistical reports.

The balance can store up to 99 accumulation records.

1. **Set the Statistic Mode in Parts Counting Mode.**
 - Statistics can be determined Manual mode or automatically.
 - **Automatic Mode:** Weights are automatically recorded when stable.
 - To use the statistics function, enable it first in the Item Settings.
 - Add items and press the **+** Statistic button to accumulate the statistical data.
 - Tap the statistics number to view the statistical results. (See Figures 30-31)

Figure 30



Figure 31

2. **Statistics Report View:**
 - Press the button displaying the number of parts counted (button to the left of +) to view statistical results. The report content includes Start Time, Number of Samples, Total, Average, Maximum, Minimum, Range, Standard Deviation, Relative Deviation % and individual item weighing values. (See Figures 32-33)

Figure 32



Figure 33

3. **Print Report**

After reviewing the report, users can choose from several actions:

- Delete the statistics and restart the percent weighing process.
- Print the report to a printer or PC based on the Print Settings.
- Export the report as a PDF file to a USB flash drive.

LIBRARY AND PRINTING SETTINGS

NOTE: For details, see Library on page 45 and Printing Settings on page 46.

CHECK COUNTING

Checks if the current sample pieces are within tolerances (e.g. an over and under limit).

NOTE: Set the sample size and average piece weight (APW) before starting check counting application.

1. Select "Check Counting" in the lower left side of the display.
2. Press ***TARE*** or ***ZERO*** to start weighing.
3. Place objects on the pan to check if the number of pieces are within the tolerances. When stable, the * appears. The default over and under limit pieces is displayed.

REFERENCE FIELD

- **Gross:** Displays the gross weight when the container is tared.
- **Net:** Displays the net weight when the container is tared.
- **Tare:** Displays the tare value when the container is tared.
- **APW:** Displays the average piece weight.
- **Over Limit:** Displays set maximum limit.
- **Under Limit:** Displays set minimum limit.
- **Check Indicator:** Displays the proportion of the current load relative to the total capacity. Each segment corresponds to 5% of total capacity. (See Figure 34)
 - **Green:** Indicates within set limit.
 - **Red:** Indicates over set limit.
 - **Orange:** Indicates under set limit.

Figure 34

- **Status Indicator:** Displays the color-coded status corresponding with the sample weight. (See Figure 35)

Figure 35

APPLICATION BUTTONS

BUTTON	FUNCTION
SAMPLE SIZE	Input sample size.
SET APW	Place the sample on the pan. Balance will calculate average piece weight (APW).
INPUT APW	Input an established APW.
0	Updates to display the number of parts counted. Press to review the statistical results.
+	Press to manually accumulate weight. In automatic statistics mode, this button will be inactive.

ITEM SETTINGS

Press the ***ITEM SETTINGS*** button to change the application settings.

Auto Optimization

Automatically optimize APW to improve counting accuracy.

Statistics

Choose statistics mode.

- **Automatic:** Weights are automatically recorded when stable.
- **Manual:** User must press the **+** button to record statistical data.

Audible Signal

The balance beeper will sound to alert the user of the check status.

Available Settings: Off, Under, Accept, Over, Under and Over

Secondary Weighing Unit

Choose between pre-set weighing units and two customized units.

Auto Tare

Choose whether to automatically tare the container value. Default setting is off.


When auto tare is enabled, the display shows the message "Place container on the pan". After a container is placed on the weighing pan, its weight will be stored in the balance. The net weight and gross weight values will then be displayed in the reference field.

SET A KNOWN AVERAGE PIECE WEIGHT

Each time a new type of part is counted, the nominal weight of one piece must be established using a small quantity of pieces.


1. Press ***INPUT APW*** and enter the desired APW value. (See Figure 36)

Figure 36

2. Press ***ENTER*** to save the value. (See Figure 37)

Figure 37

SET A DERIVED AVERAGE PIECE WEIGHT

1. Press ***SAMPLE SIZE*** and enter the desired sample size. The sample size can be 1-10,000 pieces. The default sample size is 10.
2. Press ***ENTER*** to save the value.
3. Place reference weight on the pan and press ***ACCEPT*** to establish a new APW. The home screen will display the new APW. (See Figures 38-39)

NOTE: Once a sample size is changed, the balance will immediately recalculate APW screen.

Figure 38




Figure 39

ESTABLISH OVER LIMIT AND UNDER LIMIT AND CHECK THE RESULT

1. Over limit and under limit are set in library. The main screen will display the check status of the current sample in the library.
2. If the check status is either under or over, the lights will turn red to indicate this.
3. When the piece value exceeds the over limit, the screen will display "Over". When the piece value is below the under limit, the screen will display "Under". When the piece value is in the accept range, the screen will display "Accept". Example: The home screen shows 30 pieces in the accept range. (See Figure 40)


Figure 40

ACCUMULATE DATA

1. Set statistics mode in item settings. (See Figure 41)

Figure 41

2. Press the button displaying the number of parts counted (button to the left of **+**) to view statistical results. The report content includes Start Time, Number of Samples, Total, Average, Maximum, Minimum, Range, Standard Deviation, Relative Deviation % and individual item weighing values. (See Figures 42-43)

Figure 42



Figure 43

3. After reviewing the report, users can choose from several actions:
 - Delete the statistics and restart the percent weighing process.
 - Print the report to a printer or PC based on the print settings.
 - Export the report as a PDF file to a USB flash drive.

LIBRARY AND PRINTING SETTINGS

NOTE: For details, see Library on page 45 and Printing Settings on page 46.

PERCENT WEIGHING

Displays the current weight as a percentage of a reference weight.

NOTE: Set the reference weight, reference adjust or the set reference weight before starting percent weighing application.

1. Select "Percent Weighing" in the lower left side of the display.
2. Press ***TARE*** or ***ZERO*** to start weighing.
3. Place objects on the pan. The difference between the sample and reference weight is displayed as both a weight and a percentage.
4. When stable, the * appears. The balance will display the percentage of the reference weight. The reference weight value, the set of reference factor and difference factor will display in the reference field.

REFERENCE FIELD

- **Gross:** Displays the gross weight when the container is tared.
- **Net:** Displays the net weight when the container is tared.
- **Tare:** Displays the tare value when the container is tared.
- **Ref. Weight:** Displays the set reference weight.
- **Ref. Adjust**: Displays the set value adjusting the reference weight.
- **Difference Percent:** Displays the percent difference between the sample weight and the predetermined reference weight.

APPLICATION BUTTONS

BUTTON	FUNCTION
SET REF. WEIGHT	Place the sample on the pan. The balance will store the reference weight.
REFERENCE WEIGHT	Input the reference weight.
REFERENCE ADJUST	Set a known value used to adjust the reference weight.
0	Updates to display the number of percent weighing actions. Press to review the statistical results.
+	Press to manually accumulate weight. In automatic statistics mode, this button will be inactive.

ITEM SETTINGS

Press the ***ITEM SETTINGS*** button to change the application settings.

Statistics

Choose statistics mode.

- **Automatic:** Weights are automatically recorded when stable.
- **Manual:** User must press the **+** button to record statistical data.

Secondary Weighing Unit

Choose between pre-set weighing units and two customized units.

Auto Tare


Choose whether to automatically tare the container value. Default setting is off.

When auto tare is enabled, the display shows the message "Place container on the pan". After a container is placed on the weighing pan, its weight will be stored in the balance. The net weight and gross weight values will then be displayed in the reference field.

SET A KNOWN REFERENCE WEIGHT

1. Press ***REFERENCE WEIGHT*** and enter the desired reference weight. The reference weight can be set within the range of 100d to the maximum capacity of the balance.
2. Press ***ENTER*** to save the value. The main screen will display the new reference weight. (See Figure 44)


Figure 44

SET A DERIVED REFERENCE WEIGHT

1. Place reference weight on the pan, and press ***ACCEPT*** to establish a new reference weight.
2. The main screen will display the new reference weight. (See Figure 45)

Figure 45

SET A REFERENCE ADJUST


1. Press the ***REFERENCE ADJUST*** button and set the adjust percentage. Adjust percentage can be 0.1-440%. (See Figure 46)

Figure 46

2. Press ***ENTER*** to save the value. If a library is activated, the new reference weight will overwrite the existing value in the library. The main screen will display the reference adjust factor. (See Figure 47)


Figure 47

ACCUMULATE DATA

1. Set statistics mode in item settings. (See Figure 48)

Figure 48

2. Press the button displaying the number of percent weighing actions (button to the left of **+**) to view statistical results. The report content includes Start Time, Number of Samples, Total, Average, Maximum, Minimum, Range, Standard Deviation, Relative Deviation % and individual item weighing values. (See Figures 49-50)

Figure 49



Figure 50

3. After reviewing the report, users can choose from several actions:
 - Delete the statistics and restart the percent weighing process.
 - Print the report to a printer or PC based on the print settings.
 - Export the report as a PDF file to a USB flash drive.

LIBRARY AND PRINTING SETTINGS

NOTE: For details, see Library on page 45 and Printing Settings on page 46.

CHECK WEIGHING

Checks if the current weight of sample pieces is within tolerances (e.g. an over and under limit, nominal weight and +/– tolerances).

NOTE: Set the over and under limits and/or nominal weight and tolerances before starting check weighing application.

1. Select "Check Weighing" in the lower left side of the display.
2. Press ***TARE*** or ***ZERO*** to start weighing.
3. Place objects on the pan to check if the number of pieces is within the tolerances. When stable, the * appears. The default over and under limit pieces is displayed.

REFERENCE FIELD

- **Gross:** Displays the gross weight when the container is tared.
- **Net:** Displays the net weight when the container is tared.
- **Tare:** Displays the tare value when the container is tared.
- **Over Limit:** Displays set maximum limit.
- **Under Limit:** Displays set minimum limit.

- **Check Indicator:** Displays the proportion of the current load relative to the total capacity. Each segment corresponds to 5% of total capacity. (See Figure 51)
 - **Green:** Indicates within set limit.
 - **Red:** Indicates over set limit.
 - **Orange:** Indicates under set limit.

Figure 51

- **Status Indicator:** Displays the color-coded status corresponding with the sample weight. (See Figure 52)

Figure 52

ITEM SETTINGS

Press the ***ITEM SETTINGS*** button to change the application settings.


Limits Setting Mode

Check result with different limits settings.

NOTE: Limits are set in the library. The main screen will display the check status of the current sample in the library.


1. In ***ITEM SETTINGS***, select limits setting mode: Under and Over; Nominal, wt. Tolerance; Nominal, % Tolerance. (See Figure 53)

Figure 53


2. In ***SAMPLE LIBRARY***, enter the required data:
 a. **Over and Under:** Set over and under limits. (See Figure 54)

Figure 54


 b. **Nominal, wt. Tolerance:** Set nominal weight, + tolerance and – tolerance. (See Figure 55)

Figure 55

 c. **Nominal, % Tolerance:** Set nominal weight, +% tolerance and –% tolerance. (See Figure 56)

Figure 56

Display

Choose how result will display on the main screen.

- **Weight:** The actual weight will display on the main screen with the under/accept/over indicator above.
- **Check Status:** "Accept", "Under", or "Overload" will display in large text on the main screen along with the under/accept/over indicator. (See Figures 57-59)

Figure 57


Figure 58



Figure 59

Audible Signal

The balance beeper will sound to alert the user of the check status.

Available Settings: Off, Under, Accept, Over, Under and Over.

Primary Weighing Unit

Choose weighing unit to display on the main screen. Default setting is grams.

Secondary Weighing Unit

Choose between pre-set weighing units and two customized units.

LIBRARY AND PRINTING SETTINGS

NOTE: For details, see Library on page 45 and Printing Settings on page 46.

DYNAMIC WEIGHING

Set to weigh objects that are not stable, such as animals.

NOTE: Set the start mode, weighing average time and auto print result before starting dynamic weighing application

1. Select "Dynamic Weighing" in the lower left side of the display.
2. Press ***TARE*** or ***ZERO*** to start weighing.
3. Place the moving sample on the pan. The balance will immediately calculate the average weight.
4. The balance will start a countdown based on the pre-set time configuration. The moving mouse icon will stop at the end of the process and the final weight will be displayed in orange.

REFERENCE FIELD

- **Gross:** Displays the gross weight when the container is tared.
- **Net:** Displays the net weight when the container is tared.
- **Tare:** Displays the tare value when the container is tared.
- **Sample Name:** Displays the set sample name.
- **Sample ID:** Displays the set sample ID.
- **Time:** Displays set averaging time.

APPLICATION BUTTON

BUTTON	FUNCTION
START	Begin dynamic weighing.

ITEM SETTINGS

Press the ***ITEM SETTINGS*** button to change the application settings.

Start Mode

- **Automatic:** User does not have to press ***START*** after weighing second sample.
- **Manual:** User must press ***START*** for each sample weighing.

Auto Print Result

Averaged result will be printed immediately. User will not need to press ***PRINT*** on main screen.

Primary Weighing Unit

Choose weighing unit to display on the main screen. Default setting is grams.

Auto Tare

Choose whether to automatically tare the container value. Default setting is off.

When auto tare is enabled, the display shows the message "Place container on the pan". After a container is placed on the weighing pan, its weight will be stored in the balance. The net weight and gross weight values will then be displayed in the reference field.

SET WEIGHING AVERAGING TIME

NOTE: Weighing averaging time is set in the library. The main screen will display the averaging weight of the current sample.

- In ***SAMPLE LIBRARY***, enter the desired averaging time. The averaging time can be 3-120 seconds.

LIBRARY AND PRINTING SETTINGS

NOTE: For details, see Library on page 45 and Printing Settings on page 46.

TOTALIZATION

Used to sum the weights of multiple samples and report the statistical data for the series of samples.

NOTE: Set the start mode before starting totalization application.

1. Select "Totalization" in the lower left side of display.
2. Press ***TARE*** or ***ZERO*** to start weighing.
3. Place each sample on the pan. When stable, the * appears. Press ***ACCUMULATE*** to add up the weights. The total weight of the samples will be displayed in the reference field. (See Figures 60-61)

Figure 60



Figure 61

REFERENCE FIELD

- **Gross:** Displays the gross weight when the container is tared.
- **Net:** Displays the net weight when the container is tared.
- **Tare:** Displays the tare value when the container is tared.
- **Total:** Displays the total weight of all samples weighed.
- **Samples:** Displays the number of samples weighed.
- **Average:** Displays the average weight of all samples.


APPLICATION BUTTONS

BUTTON	FUNCTION
ACCUMULATE	Add up the sample weights.
RESULT	Review accumulated data.

ITEM SETTINGS

Press the ***ITEM SETTINGS*** button to change the application settings. (See Figure 62)

Figure 62

Start Mode

- **Automatic:** User does not have to press ***ACCUMULATE*** after weighing second sample.
- **Manual:** User must press ***ACCUMULATE*** for each sample weighing.

Primary Weighing Unit

Choose weighing unit to display on the main screen. Default setting is grams.

Auto Tare

Choose whether to automatically tare the container value. Default setting is off.

When auto tare is enabled, the display shows the message "Place container on the pan". After a container is placed on the weighing pan, its weight will be stored in the balance. The net weight and gross weight values will then be displayed in the reference field.

ACCUMULATE DATA

Balance can store up to 99 accumulation records.

1. Press ***RESULT*** to view all data. The report content includes Start Time, Number of Samples, Total, Average, Maximum, Minimum, Range, Standard Deviation, Relative Deviation % and individual item weighing values. (See Figures 63-64)

Figure 63



Figure 64

2. After reviewing the report, users can choose from several actions:
 - Delete the statistics and restart the percent weighing process.
 - Print the report to a printer or PC based on the print settings.
 - Export the report as a PDF file to a USB flash drive.

LIBRARY AND PRINTING SETTINGS

NOTE: For details, see Library on page 45 and Printing Settings on page 46.

FORMULATION

Used to combine various elements in proportionate amounts.

Both free formulate and recipe based modes allow user to increase the volume of the formulation, enabling the dosage form to achieve the desired size and weight. This is particularly important for formulations containing small amounts of active pharmaceutical ingredients (APIs).

NOTE: Set the filler weight before starting formulation application.

1. Select "Formulation" in the lower left side of the display.
2. Press ***TARE*** or ***ZERO*** to start weighing.
3. Place the ingredients on the pan. The total weight of the samples will be displayed in the reference field.

REFERENCE FIELD

- **Gross:** Displays the gross weight when the container is tared.
- **Net:** Displays the net weight when the container is tared.
- **Tare:** Displays the tare value when the container is tared.
- **Total Weight:** Displays the total weight of the formulation (free formulation only).
- **Diff (wt.):** Displays weight difference between actual weight of ingredient and the ingredient weight listed in the recipe (recipe based only).
- **Diff (%):** Displays percentage difference between actual weight of ingredient and the ingredient weight listed in the recipe (recipe based only).
- **Target:** Displays the target weight based on the recipe library (recipe based only).

APPLICATION BUTTONS

BUTTON	FUNCTION
INPUT FILLER	Input filler value.
SET FILLER	Place filler on the pan and press to accept the value.
CREATE RECIPE (Recipe Based Only)	Create a new recipe.
VIEW RESULT	View the formulation results once the process has finished.
START	Start formulation.

ITEM SETTINGS


Press the ***ITEM SETTINGS*** button to change the application settings.

Formulation Mode

Choose method of formulation. (See Figure 65)

- **Free Formulate:** User freely adds ingredients.
- **Recipe Based:** User must enter the formula first, followed by individual ingredients.

Figure 65

Filler

Set filler on or off.

Primary Weighing Unit

Choose weighing unit to display on the main screen. Default setting is grams.

Secondary Weighing Unit

Choose between pre-set weighing units and two customized units.

Auto Tare


Choose whether to automatically tare the container value. Default setting is off.

When auto tare is enabled, the display shows the message "Place container on the pan". After a container is placed on the weighing pan, its weight will be stored in the balance. The net weight and gross weight values will then be displayed in the reference field.

FREE FORMULATE

1. Press ***ITEM SETTINGS***. In "Formulation Mode", select "Free Formulate". (See Figure 66)

Figure 66

2. If "Filler" option has been enabled in settings, user must enter a filler value before creating a free recipe. User can either press ***INPUT FILLER*** button on the main screen and set the value or place the filler on the weighing pan and store the displayed value as the filler weight. (See Figures 67-68)

Figure 67




Figure 68


3. Place a sample on the weighing pan.
4. Press the ***ACCEPT*** to add the ingredient's weight to the formulation. Repeat for remaining ingredients. (See Figure 69)

Figure 69

5. User can choose to add the free formulation to recipe for future use in recipe based formulation. Once the process is completed, a result screen will appear with a button to review the data. Print the result or transfer the data to a PDF file using a USB flash drive. (See Figure 70)

Figure 70


RECIPE BASED

NOTE: The system can store up to 99 ingredients in a recipe. Up to 25 recipes can be set up.


1. Press ***ITEM SETTINGS***. In "Formulation Mode", select "Recipe Based". (See Figure 71)

Figure 71


2. Press ***CREATE RECIPE*** and create a new recipe name. (See Figure 72)

Figure 72

3. Select ***NEW*** to activate a recipe. Press ***EDIT*** to go to the next step. (See Figure 73)

Figure 73

4. If "Filler" option has been enabled in settings, user must enter a filler value before creating a recipe. User can either press ***INPUT FILLER*** button on the main screen and set the value or place the filler on the weighing pan and store the displayed value as the filler weight.
5. Press ***ADD NEW*** to create ingredient names and add their corresponding weights in grams. (See Figures 74-75)

Figure 74




Figure 75

6. Use the ***UP*** and ***DOWN*** buttons to adjust the sequence order of the ingredients.
7. Press ***SAVE*** before exiting the recipe screen. Once the recipe has been saved, the user can press the ***START*** button to begin the formulation process.
8. Weigh the ingredients in order specified in the active recipe. Place the sample on the weighing pan and press ***ACCEPT*** to confirm the weight. Follow the on-screen guide messages to repeat for each ingredient.
9. Once the process is completed, a result screen will appear with a button to review the data. Print the result or transfer the data to a PDF file using a USB flash drive. (See Figure 76)

Figure 76

LIBRARY AND PRINTING SETTINGS

NOTE: For details, see Library on page 45 and Printing Settings on page 46.

DIFFERENTIAL

Used to calculate the difference in weights of multiple samples taken at different times.

NOTE: Set the sample groups number before starting differential weighing application. The initial weight, final weight and difference weight will be displayed in the reference field.

1. Select "Differential" in the lower left side of the display.
2. Press ***TARE*** or ***ZERO*** to start weighing.
3. Input desired sample group number. Sample group number can be 1-99.
4. Press ***START*** and place the sample on the pan. Press ***ACCEPT*** to confirm the initial weight.
5. Follow the on-screen messages, repeating step 4 until initial weights of all samples have been captured.
6. Weigh the samples again in the same order to capture the final weights.
7. When complete, the display will show the results. Print the result or transfer the data to a PDF file using a USB flash drive.

REFERENCE FIELD

- **Gross:** Displays the gross weight when the container is tared.
- **Net:** Displays the net weight when the container is tared.
- **Tare:** Displays the tare value when the container is tared.
- **Initial Weight:** Displays the initial weight of the sample.
- **Final Weight:** Displays the final weight of the sample.
- **Diff (wt.):** Displays the difference between the initial and final weights.

APPLICATION BUTTONS

BUTTON	FUNCTION
SAMPLE GROUPS	Set the number of sample groups requiring both initial and final weight measurements.
VIEW RESULT	View the differential results once the process has finished.
START	Start differential weighing.
ACCEPT	Accept the weighing on the pan.
END PROCESS	Complete the differential weighing process.


ITEM SETTINGS

Press the ***ITEM SETTINGS*** button to change the application settings.

Sample Sequence

- **Off:** Weigh the initial weight of all samples, then weigh the final weight of all the samples. Example: Weigh samples 1, 2, 3, 4 and 5 for their initial weights. Then weigh samples 1, 2, 3, 4 and 5 again for their final weights.
- **A set of the initial and final weight:** Weigh the initial and final weights of each sample as a pair. Example: Weigh sample 1 for its initial weight. Then weigh sample 1 again for its final weight. Repeat for all remaining samples. (See Figure 77)

Figure 77

Reciprocal Proportion

Choose how difference weight will be calculated.

- **Off:** Difference weight = final weight – initial weight
- **On:** Difference weight = initial weight – final weight

Absolute Value

The difference result will be positive even if initial weight is heavier than final weight.

On: Difference weight = | Final weight – initial weight |


PEAK HOLD

Captures the maximum weight in a series of weighing.

NOTE: Select the hold mode before starting peak hold weighing application.

1. Select "Peak Hold" in the lower left side of display.
2. Press ***TARE*** or ***ZERO*** to start weighing. (See Figure 78)

Figure 78

REFERENCE FIELD

- **Gross:** Displays the gross weight when the container is tared.
- **Net:** Displays the net weight when the container is tared.
- **Tare:** Displays the tare value when the container is tared.
- **Mode:** Displays selected hold mode.
- **Start Mode:** Displays selected start mode.
- **Sample Name:** Displays selected sample name.

APPLICATION BUTTONS

BUTTON	FUNCTION
START (Peak Hold Only)	Begin the peak hold weighing.
STOP (Peak Hold Only)	Reset the peak hold.
HOLD (Display Hold Only)	Record the last weight.
CANCEL (Display Hold Only)	Reset the display hold.

ITEM SETTINGS

Press the ***ITEM SETTINGS*** button to change the application settings.

Mode

Choose hold mode. Default setting is peak hold.

- **Peak Hold:** Balance result screen will hold the highest weighing value.
- **Display Hold:** Balance result screen will hold the last weighing value.

Start Mode

Choose start mode. Default setting is manual.

- **Manual:** User must press a button to start weighing, hold a weight value and reset the held value.
- **Semi-Automatic:** Balance begins weighing when an object is placed on the pan. User must reset the held value manually.
- **Automatic:** Balance begins weighing when an object is placed on the pan, and held value is automatically reset when the object is removed.

Stable Weight Only

Choose whether to capture weight only when weight is stable. Default setting is off.

Primary Weighing Unit

Choose weighing unit to display on the main screen. Default setting is grams.

Secondary Weighing Unit

Choose between pre-set weighing units and two customized units.

Auto Tare


Choose whether to automatically tare the container value. Default setting is off.

When auto tare is enabled, the display shows the message "Place container on the pan". After a container is placed on the weighing pan, its weight will be stored in the balance. The net weight and gross weight values will then be displayed in the reference field.

PEAK HOLD

1. Press ***ITEM SETTINGS***. In "Mode", select "Peak Mode".
2. Select start mode (Manual, Semi-Automatic, Automatic). (See Figure 79)

Figure 79

3. Place a sample on the weighing pan. Press ***START*** to check the weight. Press ***STOP*** to restart the peak hold mode. (See Figures 80-81)

Figure 80




Figure 81

DISPLAY HOLD


1. Press ***ITEM SETTINGS***. In "Mode", select "Display Hold". (See Figure 82)

Figure 82

2. Place a sample on the weighing pan. Press ***HOLD*** to record the last weight. Press ***CANCEL*** to restart the display hold mode. (See Figure 83)

Figure 83

LIBRARY AND PRINTING SETTINGS

NOTE: For details, see Library on page 45 and Printing Settings on page 46.

DENSITY DETERMINATION


Used to determine the density of a solid or a liquid.

NOTE: A density determination kit (not included) must be installed on the balance. For installation instructions, refer to the density kit manual. The balance software includes a built-in reference density table for water at temperatures ranging from 50°F to 86°F. Before attempting density measurements, please review this entire section thoroughly

NOTE: Refer to *USER HELP* before setting up the process.

1. Select "Density Determination" in the lower left side of the display.
2. Press ***TARE*** or ***ZERO*** to start weighing. (See Figure 84)

Figure 84

REFERENCE FIELD

- **Liquid Density:** Displays the current auxiliary liquid density value. Default setting is according to the distilled water temperature.
- **Water Temp.:** Displays the auxiliary liquid temperature. Default setting is according to the distilled water temperature. The water temperature range is 50° to 86°F.
- **Oil Density:** Displays the density of the selected oil (porous material only).
- **Dry Weight:** Displays the weight of the dry sample (porous material only).
- **Oiled Weight:** Displays the weight of the oiled sample (porous material only).
- **Liquid Weight:** Displays the weight of the sample in liquid.
- **Air Weight:** Displays the sample weight in the air (solid and liquid materials only).
- **Sinker Volume:** Displays sinker volume (liquid material only).

APPLICATION BUTTONS

BUTTON	FUNCTION
USER HELP	Outlines the procedures for determining the density of solid materials.
START	Begin density determination process.
CANCEL	Cancel the current process.
ACCEPT	Record the current weight on the display.

ITEM SETTINGS

Press the ***ITEM SETTINGS*** button to change the application settings.

g/cm³

Choose the resolution of the density value.

- **H-12947 (EXR224):** 0.1g/cm³, 0.01g/cm³, 0.001g/cm³, 0.0001g/cm³, 0.00001g/cm³. Default setting is 0.001g/cm³.
- **H-12948 (EXR4202) and H-12949 (EXR6202):** 0.1g/cm³, 0.01g/cm³. Default setting is 0.01g/cm³.

Density Type

Set the sample's density type.

Available Settings: Solid, Liquid. Default setting is solid.

Liquid Type

Set the liquid type.

Available Settings: Water, Other Liquid. Default setting is water.

Porous Materials

Set whether material is porous.

Available Settings: On, Off. Default setting is off.

Auto Print Result

Density result will be printed immediately. User will not need to press ***PRINT*** on main screen.

Auto Sample

Auto process one sample after another.

Auto Tare

Choose whether to automatically tare the container value. Default setting is off.

When auto tare is enabled, the display shows the message "Place container on the pan". After a container is placed on the weighing pan, its weight will be stored in the balance. The net weight and gross weight values will then be displayed in the reference field.

DENSITY DETERMINATION FOR SOLID MATERIAL

The principle of measuring the density of a solid using a balance is based on the definition of density and the measurement of mass and volume. Density is defined as the mass of an object divided by its volume.


When measuring the density of a solid, distilled water is generally used as the medium. This is because the density of distilled water is known (it reaches its maximum density of 1 g/cm³ at 39.2°F), and it does not chemically react with most solid materials.

1. Set up the density kit on the balance. Follow the instructions included with kit.
2. In ***ITEM SETTINGS***, confirm the following are selected:
 a. **Density Type:** Solid
 b. **Liquid Type:** Water

NOTE: If using another liquid, change the liquid density in sample library.


 c. **Porous Material:** Off
3. In ***SAMPLE LIBRARY***, set water temperature. Default setting is 68°F. (See Figure 85)

Figure 85


4. Press ***BACK*** to return to the main screen.
5. Prepare the sample. Ensure the solid material sample is clean and dry.
6. Press ***START*** to weigh the dry material sample in air. Balance will record its mass. (See Figure 86)

Figure 86


7. Push the sample down into the liquid until it is fully submerged. Ensure water fills all cavities and no bubbles are trapped.
8. After the sample is fully saturated, weigh the sample again. Press ***ACCEPT***. (See Figure 87)

Figure 87

9. Density result will be shown on the display. Print the result or transfer the data to a PDF file using a USB flash drive. (See Figure 88)


Figure 88

DENSITY DETERMINATION FOR POROUS MATERIAL


1. Set up the density kit on the balance. Follow the instructions included with kit.
2. In ***ITEM SETTINGS***, confirm the following are selected:
 a. **Density Type:** Solid
 b. **Liquid Type:** Oil
 c. **Porous Material:** On
3. In ***SAMPLE LIBRARY***, set the oil density. Default setting is 0.8000 g/cm³. (See Figure 89)

Figure 89

4. Press ***BACK*** to return to the main screen.
5. Select an appropriate oil that has good wettability with the porous material. Commonly used oils include kerosene, light machine oil, etc.
6. Prepare the sample. Ensure the porous material sample is clean and dry.
7. Press ***START*** to weigh the dry porous material sample in air. Balance will record its mass. (See Figure 90)


Figure 90

8. Push the sample down into the oil until it is fully submerged. Ensure oil fills all open pores in the material.
9. After the sample is fully saturated, weigh the oiled sample in air. Press ***ACCEPT***. (See Figure 91)


NOTE: The weight of the oiled sample must be less than the weight of the sample in air.

Figure 91


10. After the sample is fully saturated, weigh the sample again while still submerged in the oil. Press ***ACCEPT***. (See Figure 92)

Figure 92

11. Density result will be shown on the display. Print the result or transfer the data to a PDF file using a USB flash drive. (See Figure 93)

Figure 93


DENSITY DETERMINATION FOR LIQUID MATERIAL

1. Set up the density kit on the balance. Follow the instructions included with kit.

NOTE: A sinker is required, as it serves as a standard reference for determining liquid density.


2. In ***SAMPLE LIBRARY***, set sinker volume. Default setting is 10.0 mL. (See Figure 94)

Figure 94


3. Press ***BACK*** to return to the main screen.
4. Prepare the sample. Ensure the liquid is free of bubbles or gas.
5. Press ***START*** to weigh the sinker in air. Balance will record its mass. Press ***ACCEPT***. (See Figure 95)

Figure 95


6. Submerge the sinker completely in the liquid. Press ***ACCEPT***. (See Figure 96)

Figure 96

7. Density result will be shown on the display. Print the result or transfer the data to a PDF file using a USB flash drive. (See Figure 97)

Figure 97

LIBRARY AND PRINTING SETTINGS

NOTE: For details, see Library on page 45 and Printing Settings on page 46.

WEIGH BELOW

The balance is equipped with a weigh below hook for weighing below the balance.


NOTE: Before turning balance over, remove the weighing pan, fixing ring, EMC plate and draft shield (if present) to prevent damage.

CAUTION! Do not place the balance on the pan support cone or load cell pins.

To use this feature, remove power from the balance, then remove the protective cover for the weigh below opening (two screws). The protective cover is reversible for easy storage.

The balance can be supported using lab jacks or any other convenient method. Ensure the balance is level and secure. Power on the balance, then use a string or wire to attach items to be weighed. (See Figure 98)

Figure 98

LIBRARY

- In library menu, users can review a list of generated libraries. Libraries can be imported and exported as needed.
- Total memory library is displayed in the upper right corner. Example: 9/1% means that 9 libraries have been created and occupy approximately 1% of total library capacity. When memory usage reaches 80%, balance will display a message prompting the user to delete or back up library records.
- **Total Capacity: Approximately 1,000 library records.**


LIBRARY MENU

APPLICATION MODE	LIBRARY FUNCTION	DATA SPECIFIC TO APPLICATION MODE	SAMPLE DATA
Basic Weighing	Yes	Minimum Weight	Sample Name Sample ID Batch ID Lot ID Project Name 10 x Customer ID
Parts Counting	Yes	APW, Sample Size	
Check Counting	Yes	APW, Sample Size, Over Limit, Under Limit	
Percent Weighing	Yes	Reference Weight, Reference Adjust	
Check Weighing	Yes	Over Limit, Under Limit, Nominal Weight, + Weight Tolerance, – Weight Tolerance, + % Tolerance, – % Tolerance	
Dynamic Weighing	Yes	Average Time	
Totalization	Yes	----------	
Peak Hold	Yes	----------	
Density Determination	Yes	Sinker Volume, Water Temperature, Liquid Density, Oil Density	
Formulation	Recipe	Recipe Name, Item Name, Item Weight	----------
Differential	No Library		

CREATE A LIBRARY


1. Press ***+ NEW LIBRARY*** to create a new library. (See Figure 99)
2. When a new library is created, it remains inactive. User must activate the library by checking the box.
3. After a library item is activated, the main screen will show the activated library number.

Figure 99

4. Other sample information can be configured in the library, including sample name, sample ID, batch ID, lot ID, project ID and the creation of 10 more customized IDs.
5. Press the trash can icon to delete a library. (See Figure 100)

Figure 100

IMPORT LIBRARY

1. Save a "library.db" file to a USB flash drive. Connect the USB drive to the balance. Once connected, the ***IMPORT*** button will become active.
2. Press ***IMPORT*** to initiate the import process.

NOTE: Library will be renamed if the library being imported has the same name as an existing library.

EXPORT LIBRARY

1. Connect the USB flash drive to the balance. Once connected, the ***EXPORT*** button will become active.
2. Press the ***EXPORT*** button to download the "library.db" file to the USB flash drive.

PRINTING SETTINGS

- User can configure printing settings based on specific application requirements.

NOTE: There is 40-second timeout control for print while stable requirement. If the unstable condition continues over 40 seconds, balance will return to previous display.

PRINT CONTENT

Specify the content of the printed data. Enable or disable content as desired. (See Figures 101-102)

Figure 101



Figure 102

Available Print Contents: Print Review, Select All, Date and Time, Balance ID, Balance Name, User ID, First Name/Last Name, Project Name, Application Mode, Sample Name, Sample ID, Batch ID, Lot ID, Customized IDs, Result, Gross, Net, Tare, Signature Line, Verified Line, Feeds Lines (1, 4, 10)

CONNECT TO PRINTER

Numeric Only

- **On:** Print only the numeric weight value.
- **Off:** Print all the weighing values enabled in Print Content.

Single Header Only

When enabled, header line will only print once every 24 hours.

Header Line Contents: Date and Time, Balance ID, Balance Name, Project Name, Application Mode, Sample Name, Sample ID, Batch ID, Lot ID, Customized IDs

Manual Print

- **Stable Only:** Print only the stable value.
- **All Values:** Print all values.

Auto Print

- **Off:** Disable auto print.
- **On Stability:** Printing occurs each time stability criteria are met.
- **Interval:** Printing occurs at the defined time interval. The interval can be set from 1-3,600 seconds.
- **Continuous:** Printing occurs continuously.

When "On Stability" is selected, set the conditions for printing.

- **Load:** Prints when displayed load is stable.
- **Load Zero:** Prints when the displayed load or zero reading is stable.

CONNECT TO PC

Output Format

Choose output format based on various output strings.

- OHAUS
- SICS: MT-SICS
- ST: ST-SICS (10)

Manual Print

- **Stable Only:** Print only the stable value.
- **All Values:** Print all values.

Auto Print

- **Off:** Disable auto print.
- **On Stability:** Printing occurs each time stability criteria are met.
- **Interval:** Printing occurs at the defined time interval. The interval can be set from 1-3,600 seconds.
- **Continuous:** Printing occurs continuously.

When "On Stability" is selected, set the conditions for printing.

- **Load:** Prints when displayed load is stable.
- **Load Zero:** Prints when the displayed load or zero reading is stable.

DATA TO EXCEL

Balance supports HID connection to a computer without the need for drivers.

Format

DD-MMM-YYYY	Weight	Units
15-May-2025	100.0000	g

Manual Print

- **Stable Only:** Print only the stable value.
- **All Values:** Print all values.

Auto Print

- **Off:** Disable auto print.
- **On Stability:** Printing occurs each time stability criteria are met.
- **Interval:** Printing occurs at the defined time interval. The interval can be set from 1-3,600 seconds.
- **Continuous:** Printing occurs continuously.

When "On Stability" is selected, set the conditions for printing.

- **Load:** Prints when displayed load is stable.
- **Load Zero:** Prints when the displayed load or zero reading is stable.

Text To


- **Column:** Print all data into a single column in Microsoft Excel.
- **Cell:** Print all data into a single cell in Microsoft Excel.

SAVE TO USB

IMPORTANT! Do not unplug the USB drive during data transfer.

- **Off:** Disable save to USB.
- **Manual:** Save to USB manually.
- **Auto:** Save to USB automatically.
 - Once user has set auto-print interval time, data will be sent to the USB flash drive at the specified intervals.
 - PDF format is not available when in auto mode.
- **Format:** Select PDF or CSV. (See Figure 103)

Figure 103

COMMUNICATION

Define external communication methods and set printing parameters. Data can be output to a printer or PC.

Communication options include two USB Type-A ports, one USB-C port, one Ethernet port, one RS-232 port and optional Bluetooth® and Wi-Fi functionalities.

HID CONNECTIONS

Balance supports HID (Human Interface Device) connection to a computer without needing drivers. Selected language is displayed on the keyboard's spacebar.

English keyboard is the default setting. Click the key to change to French characters.

RS-232

Connects balance to peripheral devices such as printers, PCs and label printers, each with different output formats. Select the appropriate device type before transferring data.

PC and Printer Connections

Choose settings before connecting to a PC or printer.

- **Baud Rate (bits per second):** 1200, 2400, 4800, 9600, 19200, 38400, 115200. Default setting is 9600.
- **Transmission (number of data bits, parity, stop bits):** 7 Even 1, 7 Odd 1, 7 None 1, 8 None 1, 7 Even 2, 7 Odd 2, 7 None 2, 8 None 2. Default setting is 8 None 1.
- **Handshake (flow control method):** None, Xon/Xoff, Hardware. Default setting is Xon/Xoff.
- **12V Output:** Choose whether to allow RS-232 port to output 12V, which supports Bluetooth® adapter.

NOTE: The purpose of the switch is to prevent damage to external devices from the 12V output.

Label Printer Connections

All templates (including default) can be edited by OHAUS Label Designer. Download software on the OHAUS website.

Five label templates can be stored in the balance. Default template is #1-100 mm * 75 mm.

- **#1-100 mm * 75 mm:** Includes all basic information (Date and Time, Sample ID, Sample Name, Batch ID, Result, Gross Weight, Tare Weight, Net Weight, User Name).
- **#2-100 mm * 10 mm:** Includes weighing result only.
- **#3, #4, #5:** Empty templates that user can edit using OHAUS Label Designer software.

INTERFACE COMMANDS

The following commands will be acknowledged by the balance.

NOTE: Commands sent to the indicator must be terminated by a carriage return-line feed (CRLF). Data output is always terminated with a carriage return-line feed (CRLF).

The balance will return "ES" for invalid commands.

INTERFACE COMMANDS

COMMAND	FUNCTION
AUF	Auto login (only works when user management function is off).
I2	Inquiry of balance data.
I3	Inquiry of balance SW version and type definition number.
I4	Inquiry of serial number.
SIR	Send weight value immediately and repeat.
IP	Immediate print of displayed weight (stable or unstable). Can also be used to stop continuous print and interval print.
P	Print displayed weight according to "Stable Only" setting in the communication menu.
CP	Continuous print.
SPx	Print on stability (x = stable time; print if stability is achieved within this time).
P x	Interval print [x = print interval (1-3600 seconds)]. IP/P ends interval print. **NOTE: The corresponding settings in the communication menu are also modified.**
Z	Zero balance (same as pressing ***ZERO***).
ZI	Zero immediately.
@	Reboot.
T	Tare balance (same as pressing ***TARE***).
TI	Tare Immediately.
M x	Set current application mode to x (x depends on application, use application list).
U x	Set balance to unit x. Units vary by model.

INTERFACE COMMANDS CONTINUED

COMMAND	FUNCTION
ON	Brings out of standby.
OFF	Goes to standby.
SIU	Send weight value with currently displayed unit immediately.
C3	Begin internal calibration (same as trigger from calibration menu).
PSN	Print serial number.
PV	Print terminal software version and base software version.
# x Unit	Set counting APW (x) in unit. **NOTE: Must have APW stored; unit can be any unit.**
% x Unit	Set percent application reference weight (x) in unit. **NOTE: Must have reference weight stored; unit can be any unit.**
CO x Unit	Set check weighing over limit in x unit.
CU x Unit	Set check weighing under limit in x unit.
TIM	Print current time.
DAT	Print current date.
TIM x	Set time. x format: hh mm ss
DAT x	Set date. x format: mm dd yyyy
\EscP	Print weight immediately.
\EscT	Tare balance.
\EscU	Tare balance.
\EscV	Zero balance.
\EscW	External calibration.
\EscZ	Internal calibration.
\Escx1_#_	Print model.
\Escx3_#_	Print software version.

USB

Set up each USB port to a peripheral device before transferring data. Use USB-C connection to connect balance to a PC or transfer data directly to Microsoft Excel. Format can be in a column or cell. See Printing Settings on page 46.

ETHERNET

Connect to a Local Area Network (LAN). After balance has been connected, user can set up ethernet options in the communication menu.

- **DHCP:** When enabled, automatically obtains IP address when connecting to a LAN or a wireless network with an unknown IP address. Other Ethernet settings will be locked.
- **Static IP:** Disable DHCP and manually enter the IP address.

WI-FI AND BLUETOOTH®

Connect using an LMB842 USB dongle. The LM842 USB dongle complies with Bluetooth® 5.0 (Dual Mode) and IEEE 802.11ac Wi-Fi standards of operation.

Wi-Fi Connection

1. Enable "Network" function. Balance will search for available networks.
2. Enter the password to pair and connect to the Wi-Fi.
3. If the network is not visible, the user can press "Other…" to manually add the network.
 a. **DHCP:** When enabled, automatically obtains IP address when connecting to a LAN or a wireless network with an unknown IP address. Other Ethernet settings will be locked.
 b. **Static IP:** Disable DHCP and manually enter the IP address.
4. Once the balance connects to Wi-Fi, an icon will appear in the upper right corner of the main screen.

Bluetooth® Connection

1. Enable "Bluetooth®" function.
2. Enter the Bluetooth® device name to pair the balance to the device.
3. Once the balance connects to a Bluetooth® device, an icon will appear in the upper right corner of the main screen.
4. MAC address is used to identify the same Bluetooth® device name.


MAINTENANCE MENU

DIAGNOSTICS MENU

IR Sensors


1. Press "IR Sensors" to start testing. (See Figure 104)

Figure 104

2. Wave a hand over the sensor to trigger it and cause the sensor light to turn green.
3. When the sensors are working properly, the status icon on the right side will turn green with a check mark. (See Figure 105)


Figure 105

Digital Leveling Bubble

1. Press "Digital Leveling Bubble" to start testing.
2. Lift the balance to test leveling bubble movement.
3. When the function is working properly, the status icon on the right side will turn green with a checkmark. If the leveling bubble is not moving, the status icon will display a warning sign. (See Figure 106)

Figure 106


Terminal Status Light

1. Press "Terminal Status Light" to start testing.
2. The balance will cycle through all the light colors on the terminal.

SOFTWARE UPGRADE

1. Save the upgrade file to a USB flash drive.
2. Connect USB flash drive to the balance.
3. Press ***START*** to upgrade the system.
4. When upgrade is complete, the balance will reboot automatically. (See Figure 107)

Figure 107

BALANCE INFORMATION

Press ***BALANCE INFORMATION*** to display essential information about the balance. (See Figure 108)

SERVICE MENU

Service menu is locked to prevent authority changes and tampering with the balance. If users encounter issues while using the balance, contact Uline Customer Service at 1-800-295-5510.

Figure 108

SERVICE LOG FILE

LOG	# OF RECORDS	DATA RECORDED	PERMISSIONS
Failure Log	100	Number of times an overload or unload occurs.	----------
System Log	100,000	Menu changes, create/delete user, log ins/log offs, date and time changes, etc.	Only Administrator and Supervisor roles can export this data to PDFs in an un-editable format.
Calibration Log	3,000	User-performed calibrations and service-performed calibrations.	----------
Printed Data Log	100,000	Routine weighing data sent by manually pressing the print key or through a command from a PC.	Only Administrator and Supervisor roles can export this data to PDFs in an un-editable format.

FACTORY RESET

Reset any the following menus to their factory default settings.

- Reset All
- Quick Setup
- Calibration
- Balance Setup
- Application Modes
- Weighing Units
- Communication
- Library

GENERAL MAINTENANCE

CALIBRATION

Periodically verify calibration by placing an accurate weight on the balance and viewing the result. If calibration is required, perform an internal calibration.

CLEANING

WARNING! Electric shock hazard. Disconnect the equipment from the power supply before cleaning.

CAUTION! Do not use solvents, harsh chemicals, ammonia or abrasive cleaning agents.

- The housing may be cleaned with a cloth dampened with a mild detergent if necessary.
- Ensure liquid does not enter the interior of the balance.

END OF LIFE

- Before disposing of the electronic balance, ensure that all sensitive data and user information is erased from the device.
- If the balance is no longer functional, recycle it at a local electronic recycling center. Follow proper e-waste disposal regulations.
- Do not dispose of balance in regular trash. Check with local waste management authorities for guidance on properly disposing of electronic equipment.

SPECIFICATIONS

AMBIENT CONDITIONS

Location	Indoor Use Only
Altitude	6,561'
Ambient Temperature Range	50-86°F
Operating Temperature Range	41-104°F
Relative Humidity	Maximum 80% for temperatures up to 87.8°F; 50% at 104°F
Electrical Supply	12V DC, 1.5 A
Mains Supply Voltage Fluctuations	Up to +/-10% of nominal voltage
Installation Category	2
Pollution Degree	2

MATERIALS

Bottom Housing	Die-Cast Aluminum, Painted
Top Housing	Die-Cast Aluminum, Painted
Terminal	Die-Cast Aluminum
Weighing Platform	Plastic (PC)
Weighing Pan	316 Stainless Steel
In-Use Cover	Plastic (PET)
Draft Shield	Glass, Aluminum, Plastic

MODEL SPECIFICATIONS

ULINE MODEL NO.	H-12947 (EXR224)	H-12948 (EXR4202)	H-12949 (EXR6202)
Capacity	220 g	4,200 g	6,200 g
Weighing Units	g, mg, ct, N, oz, ozt, grain, dwt, mo, msg, tl H, tl S, tl T, tcl, tola, baht, Custom 1, Custom 2	g, ct, N, oz, ozt, grain, dwt, mo, msg, tl H, tl S, tl T, tcl, tola, baht, lb, kg, Custom 1, Custom 2	
Stabilization Time (Typical)	2 Seconds	≤ 1 Second	
Display	7" TFT 16.7M Color Touch Screen with Glass Overlay		
Communication	USB Host (Type A) x 2 USB Device (Type C) x 1 Ethernet (RJ45) x 1 RS-232 x 1 Optional Wi-Fi/Bluetooth® Dongle		
Power Input	12 V DC, 1.5A		
Power Supply	AC Adapter Input: 100-240 VAC, 0.5A, 50-60 Hz AC Adapter Output: 12 V DC, 1.5A		
Weighing Pan Size	Square 3.5 x 3.5"	Trapezoid 7.0 x 7.9"	
Draft Shield with Anti-Static Coating	Yes	No	
Assembled Dimensions (W x D x H)	9.1 x 15.5 x 14.6"	9.0 x 15.4 x 3.7"	
Net Weight	16.3 lbs.	10.8 lbs.	

WARNING NOTIFICATIONS

Warning notifications are shown on the upper right corner of the display.

CALIBRATION FAILURE

Balance has not calibrated successfully.

1. Press on the warning sign. The notification will appear on the screen.
2. Press the warning message and the screen will show the ***INITIATE CALIBRATION*** button.
3. Perform internal calibration or span calibration. See Calibration on page 10.

LEVELING THE BALANCE

If the digital leveling bubble is not centered, the balance needs to be leveled.

1. Press on the warning sign. The notification will appear on the screen.
2. Press the warning message and the screen will show the ***GO LEVELING*** button.
3. Before starting the motorized leveling process, remove any weight from the weighing pan.
4. Press the ***LEVELING*** button. Adjust the leveling feet according to the position of the digital leveling bubble until the bubble is centered.
5. When balance is level, digital leveling bubble will turn green.
6. If leveling system fails to level the balance correctly, move the balance to a flat surface and retry.

PASSWORD

The password has expired in accordance with the password policy.

1. Press on the warning sign. The notification will appear on the screen.
2. Press the warning message and the screen will show the ***RESET PASSWORD*** button.
3. Reset password.


NOTE: If user forgets their password, contact Uline Customer Service at 1-800-295-5510.

OPERATING ISSUES	POSSIBLE CAUSES	RECOMMENDATIONS
Balance will not turn on.	No power supply to balance.	Verify connection and voltage.
Inaccurate weighing results.	Improper calibration. Unstable environment.	Perform calibration. Move balance to a suitable location.
Cannot calibrate.	Calibration menu is locked. Unstable environment. Incorrect calibration masses.	Turn calibration menu lock off. Move balance to a suitable location. Use correct calibration masses.
Cannot change menu settings.	Sub-menu is locked.	Unlock sub-menu.
Low reference weight.	Reference weight is too small. The weight on the pan is too small to define a valid reference weight.	Increase sample size.
Invalid piece weight.	Average piece weight is too small.	Increase average piece weight.
Operation timeout.	Weight reading is not stable.	Move balance to a suitable location.
Display shows "-------".	Busy (tare, zero, printing).	Wait until completion.

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